

AUTHORIZED DISTRIBUTOR

MZ-Analysentechnik GmbH
 Barcelona-Allee 17 • D-55129 Mainz
 Tel +49 6131 880 96-0
 Fax +49 6131 880 96-20
 e-mail: info@mz-at.de
 www.mz-at.de

SEPARATION OF METHYLCLOTHIAZIDE USING NEW IMMOBILIZED COLUMN – CHIRALPAK IG

APPLICATION NOTE

INTRODUCTION

Daicel Corporation has introduced a new and unique chiral selector, CHIRALPAK® IG, to its line of Daicel immobilized chiral stationary phases (CSPs). The IG phase is an amylose polymer derivatized with phenyl moieties (chloro- and methyl-groups) in meta positions. This distinctive phase was developed from our extensive library of > 150 CSPs. Moreover, our studies have shown that the IG chiral selector has the highest selectivity of any of our immobilized CSPs.

Meta-substituted immobilized chiral selectors such as CHIRALPAK IA, IB and IC have been shown to have remarkable affinity for resolution of chiral compounds from different types of molecules. Based on laboratory results, we have found that the addition of CHIRALPAK IG forms the most effective set of four chiral columns, CHIRALPAK IA, IB, IC, and IG, to be utilized for primary screening. It is noteworthy that the four best columns are all meta-substituted.

Primary Screening Columns

- CHIRALPAK IA: Amylose *tris* (3,5-dimethylphenylcarbamate)
- CHIRALPAK IB: Cellulose *tris* (3,5-dimethylphenylcarbamate)
- CHIRALPAK IC: Cellulose *tris* (3,5-dichlorophenylcarbamate)
- CHIRALPAK IG: Amylose *tris* (3-chloro-5-methylphenylcarbamate)

EXPERIMENTAL AND DISCUSSION

A CHIRALPAK IG column, packed with 5- μ m particles, was used to develop the enantioselective separation of methyclothiazide. The mobile phase was a mixture of ammonium bicarbonate and acetonitrile. All chromatographic conditions are given in Figure 1.

SEPARATION EXAMPLE – REVERSED PHASE




Figure 1: Separation of Methyclothiazide

CHROMATOGRAPHIC CONDITIONS

Column Size:	Daicel CHIRALPAK IG 4.6 mm i.d. x 250 mm long
Mobile Phase:	20mM NH ₃ HCO ₃ aq/MeCN (pH 9 with DEA, 60/40)
Flow Rate:	1 ml/min
UV Detection:	254 nm
Column Temperature:	25° C

Note: Methyclothiazide is a thiazide diuretic that is used to treat high blood pressure (hypertension). Using CHIRALPAK IG provides fast and effective separation of methyclothiazide enantiomers



CHIRAL TECHNOLOGIES, INC.

800 North Five Points Road
 West Chester, PA 19380 USA
 Tel: 610-594-2100
 Fax: 610-594-2325
 www.chiraltech.com
 email: chiral@chiraltech.com

CHIRAL TECHNOLOGIES EUROPE

Parc d'Innovation, Bd Gonther d'Andernach
 67404 Illkirch Cedex, France
 Tel: +33 (0) 388 79 52 00
 Fax: +33 (0) 388 66 71 66
 www.chiraltech.com
 e-mail: cte@chiral.fr

