



# Analysis of Per- and Polyfluoroalkyl Substances in Aqueous Samples with CHROMABOND® WAX according to EPA Method 533

MACHEREY NAGEL application department · Dr. H. R. Wollseifen, T. Kretschmer, L. Emmerich

## Application benefits

- Successful determination of 25 Per- and Polyfluoroalkyl Substances from water samples according to EPA Method 533
- High recovery rates were achieved with a CHROMABOND® WAX SPE Column
- Fast and sensitive HPLC analysis on a NUCLEODUR® PFAS column

## MN products

### REF 7300011

CHROMABOND® WAX, 6 mL, 150 mg

### REF 760666.20

EC 100/2 NUCLEODUR® PFAS, 3 µm

### REF 760673.20

EC 50/2 NUCLEODUR® PFAS Delay

### REF 702402

Screw closure, N 9, PP, blue, center hole, silicone white/polyimide orange, 1 mm, fluorine-free

### REF 702009

Screw neck vial, N 9, 11.6 x 32.0 mm, 0.3 mL, inner cone, PP transparent

## MN application numbers

SPE: 306960

HPLC: 129380

## Keywords

EPA Method 533, PFAS, WAX, weak anion exchanger, water, LC-MS/MS, Delay column


## Introduction

In December 2019, the United States Environmental Protection Agency (US EPA) has published a method for the analysis of per- and polyfluoroalkyl substances (PFAS) [1]. The method was developed to target “short chain” PFAS (none greater than C12), including perfluorinated acids, sulfonates, fluoroteleners, and poly/perfluorinated ether carboxylic acids. This “short chain” PFAS could not be analyzed using 537.1 due to physicochemical properties. 25 PFAS should be analyzed using a solid phase extraction (SPE) and a liquid chromatography-tandem mass spectrometry (LC-MS/MS) method in drinking water. Method 533 requires SPE cartridge containing weak anion exchange, mixed-mode polymeric sorbent (polymeric backbone and a diamino ligand) and a particle size of approximately 33 µm. The SPE sorbent must have a pKa above 8 so that it remains positively charged during extraction. The use of 200 mg sorbents is recommended for the extraction of 100 mL samples.

In this application note, a solid phase extraction (SPE) using the CHROMABOND® WAX coupled with a liquid chromatography-tandem mass spectrometry (LC-MS/MS) method is presented. High recovery rates with very good reproducibility are achieved for drinking water matrices. Finally, the extracts are analyzed using HPLC-MS/MS on a NUCLEODUR® PFAS column.

| Requirements EPA Method 533  | Specifications of CHROMABOND® WAX |
|--|-----------------------------------|
| Weak anion exchange  | pass                              |
| Mixed-mode polymeric sorbent (polymeric backbone and a diamino ligand) | pass                              |
| Particle size approximately 33 µm                                      | pass                              |
| pKa above 8  | pass                              |
| Low blind level of PFAS  | pass                              |

Table 1: Matching of EPA 533 requirements with specifications of CHROMABOND® WAX.



### AUTHORIZED DISTRIBUTOR

MZ-Analysentechnik GmbH, Barcelona-Allee 17 • D-55129 Mainz  
Tel +49 6131 880 96-0, Fax +49 6131 880 96-20  
e-mail: info@mz-at.de, www.mz-at.de

# Analysis of PFAS in Aqueous Samples with CHROMABOND WAX according to EPA Method 533

## Sample pretreatment

MN Appl. No. 306960

### Solid phase extraction according to EPA 533

#### Sample preparation:

This method is applicable to aqueous samples containing up to 50 mg of suspended solids per sample. The procedure requires the preparation of the entire sample. Subsampling should be avoided whenever possible. Typical sample size is 500 mL.

1. Add ammonium acetate (1.0 g/L) to the sample. Ammonium acetate will sequester free chlorine to form chloramine.
2. Verify that the sample containing 1 g/L ammonium acetate has a pH between 6.0 and 8.0. Acetic acid may be added as needed to adjust the pH.
3. Add 20  $\mu$ L organic standard solution\* ( $\beta$  = 12.5 ng/mL in methanol for each compound) to the 250 mL water sample.

\* Contains native and isotopically labeled per- and polyfluoroalkyl substances

Column: CHROMABOND® WAX, 6 mL 150 mg (REF 7300011)

#### Conditioning:

Rinse each cartridge with 10 mL methanol, 10 mL of aqueous 0.1 M phosphate buffer. Close the valve and add 2–3 mL of phosphate buffer (pH 7.0)\* to the cartridge reservoir and fill the remaining volume with reagent water.

\* Mix 500 mL of 0.1 M dibasic sodium phosphate with approximately 275 mL of 0.1 M monobasic sodium phosphate. Verify that the solution pH is approximately 7.0

#### Sample application:

Attach the reservoir cartridges, turn on the vacuum, and begin adding the 250 mL water sample with a flow rate of 5 mL/min to the cartridge. Do not allow the cartridge to run dry before all the sample has passed through.

#### Sample bottle and cartridge rinse:

After the entire sample has passed through the cartridge, rinse the sample bottles and the transfer cartridges with aliquots of 1 g/L ammonium acetate in water and draw each aliquot through the SPE columns. Add 1 mL of methanol to the sample bottle and draw through the transfer cartridge and SPE cartridge. Draw air or nitrogen through the cartridge for 5 min at high vacuum (15–20 in. Hg).

#### Sample bottle and cartridge elution:

Rinse the sample bottles and the transfer cartridge with 5 mL of methanol with 2% ammonium hydroxide (v/v) and elute the analytes from the cartridges by pulling the 5 mL of methanol with 2% ammonium hydroxide (v/v) through the SPE column. Use a low vacuum such that the solvent exits the cartridge in a dropwise fashion. Repeat sample bottle rinse and cartridge elution with a second 5 mL aliquot of methanol with 2% ammonium hydroxide (v/v).

#### Eluent exchange:

Evaporate eluate to dryness at 40 °C under a stream of nitrogen and dissolve residue in 0.5 mL methanol.

## Analysis by HPLC-MS / MS

MN Appl. No. 129380

### Chromatographic conditions

|                  |   |
|------------------|---|
| DELAY Column     | EC 50/2 NUCLEODUR® PFAS Delay (REF 760673.20)   |
| Column           | EC 100/2 NUCLEODUR® PFAS, 3 $\mu$ m (REF 760666.20)   |
| Eluent A         | 5 mM ammonium acetate in water  |
| Eluent B         | 5 mM ammonium acetate in methanol   |
| Gradient         | hold 40% B for 1 min, in 8 min from 40% B to 95% B, hold 95% B for 3 min, in 0.1 min to 40% B, hold 40% B for 2.9 min |
| Flow rate        | 0.3 mL/min  |
| Temperature      | 40 °C   |
| Injection volume | 2 $\mu$ L   |
| MS conditions    |   |
| Acquisition mode | SRM   |
| Interface        | ESI   |
| Polarity         | negative  |
| Curtain Gas      | 30  |
| Collision Gas    | medium  |
| Ionspray Voltage | –4500 V   |
| Temperature      | 400 °C  |
| Ion Source Gas 1 | 50  |
| Ion Source Gas 2 | 60  |
| Detection Window | 60 sec  |



# Analysis of PFAS in Aqueous Samples with CHROMABOND WAX according to EPA Method 533

## MRM transitions

| Analyte   | Abbreviation | CAS number  | Q1 mass [Da] | Q3 mass [Da] | Retention time [min] |
|---|--------------|-------------|--------------|--------------|----------------------|
| Perfluoro-3-methoxypropanoic acid   | PFMPA        | 377-73-1    | 229.00       | 85.00        | 1.96                 |
| Perfluoro- <i>n</i> -butanoic acid  | PFBA         | 375-22-4    | 212.90       | 168.80       | 2.01                 |
| Perfluoro-4-methoxybutanoic acid  | PFMBA        | 863090-89-5 | 279.00       | 85.00        | 3.64                 |
| Perfluoro- <i>n</i> -pentanoic  | PFPeA        | 2706-90-3   | 262.88       | 219.00       | 3.90                 |
| Perfluoro(2-ethoxyethane)sulfonic acid  | PFEESA       | 113507-82-7 | 315.00       | 135.00       | 4.14                 |
| Perfluoro- <i>n</i> -butanesulfonic acid  | PFBS         | 375-73-5    | 298.93       | 98.90        | 4.20                 |
| Nonafluoro-3,6-dioxahexanoic acid   | NFDHA        | 151772-58-6 | 295.00       | 201.00       | 4.48                 |
| 1 <i>H</i> ,1 <i>H</i> , 2 <i>H</i> , 2 <i>H</i> -Perfluorohexane sulfonic acid                                     | 4:2FTS       | 757124-72-4 | 326.94       | 306.90       | 5.27                 |
| Perfluoro- <i>n</i> -hexanoic acid  | PFHxA        | 307-24-4    | 312.91       | 268.80       | 5.40                 |
| Perfluoropentanesulfonic acid   | PFPeS        | 2706-91-4   | 348.85       | 80.00        | 5.54                 |
| Hexafluoropropylene oxide dimer acid  | HFPO-DA      | 13252-13-6  | 284.99       | 168.70       | 5.77                 |
| Perfluoro- <i>n</i> -heptanoic acid   | PFHpA        | 375-85-9    | 362.93       | 318.80       | 6.45                 |
| Perfluoro- <i>n</i> -hexanesulfonic acid  | PFHxS        | 355-46-4    | 398.94       | 79.80        | 6.49                 |
| 4,8-Dioxa-3 <i>H</i> -perfluorononanoic acid  | ADONA        | 919005-14-4 | 376.90       | 250.70       | 6.58                 |
| 1 <i>H</i> ,1 <i>H</i> , 2 <i>H</i> , 2 <i>H</i> -Perfluorooctane sulfonic acid                                     | 6:2FTS       | 27619-97-2  | 426.93       | 406.90       | 7.24                 |
| Perfluoro- <i>n</i> -heptanesulfonic acid   | PFHpS        | 375-92-8    | 448.93       | 79.80        | 7.26                 |
| Perfluoro- <i>n</i> -octanoic acid  | PFOA         | 335-67-1    | 412.91       | 369.00       | 7.26                 |
| Perfluoro- <i>n</i> -octanesulfonic acid  | PFOS         | 1763-23-1   | 498.84       | 79.90        | 7.89                 |
| 3-Perfluoroheptyl propanoic acid  | 7:3FTCA      | 812-70-4    | 441.00       | 317.00       | 8.00                 |
| 9-Chlorohexadecafluoro-3-oxanonane-1-sulfonic acid  | 9Cl-PF3ONS   | 73606-19-6  | 530.75       | 350.70       | 8.25                 |
| Perfluoro- <i>n</i> -decanoic acid  | PFDA         | 335-76-2    | 512.84       | 468.90       | 8.49                 |
| 1 <i>H</i> ,1 <i>H</i> , 2 <i>H</i> , 2 <i>H</i> -Perfluorodecane sulfonic acid                                     | 8:2FTS       | 39108-34-4  | 526.00       | 506.80       | 8.50                 |
| Perfluoro- <i>n</i> -undecanoic acid  | PFUnDA       | 2058-94-8   | 562.80       | 518.90       | 8.95                 |
| 11-Chloroeicosafluoro-3-oxaundecane-1-sulfonic acid   | 11Cl-PF3OUdS | 763051-92-9 | 630.74       | 451.00       | 9.15                 |
| Perfluoro- <i>n</i> -dodecanoic acid  | PFDoDA       | 307-55-1    | 612.79       | 568.90       | 9.33                 |
| Surrogates  |              |             |              |              |                      |
| Perfluoro-(2,3,4- <sup>13</sup> C <sub>3</sub> )butanoic acid   | M3PFBA       |             | 216.00       | 172.00       | 1.91                 |
| Perfluoro-( <sup>13</sup> C <sub>4</sub> )butanoic acid   | M4PFBA       |             | 216.94       | 171.90       | 2.01                 |
| Sodium perfluoro-(2,3,4- <sup>13</sup> C <sub>3</sub> )butanesulfonate  | M3PFBS       |             | 301.89       | 98.90        | 4.22                 |
| Sodium 1 <i>H</i> ,1 <i>H</i> ,2 <i>H</i> ,2 <i>H</i> -perfluoro(1,2- <sup>13</sup> C <sub>2</sub> )hexanesulfonate | M2-4:2FTS    |             | 328.97       | 81.00        | 5.26                 |
| Tetrafluoro-2-heptafluoropropoxy- <sup>13</sup> C <sub>3</sub> -propanoic acid                                      | M3HPFO-DA    |             | 287.00       | 169.00       | 5.44                 |
| Perfluoro-(1,2,3,4- <sup>13</sup> C <sub>4</sub> )heptanoic acid  | M4PFHpA      |             | 366.95       | 321.80       | 6.45                 |
| Sodium perfluoro-(1,2,3- <sup>13</sup> C <sub>3</sub> )hexanesulfonate  | M3PFHxS      |             | 401.90       | 79.90        | 6.50                 |
| Perfluoro-(1,2,3,4- <sup>13</sup> C <sub>4</sub> )octanoic acid   | MPFOA        |             | 417.00       | 372.00       | 7.08                 |
| Sodium 1 <i>H</i> ,1 <i>H</i> ,2 <i>H</i> ,2 <i>H</i> -perfluoro(1,2- <sup>13</sup> C <sub>2</sub> )octanesulfonate | M2-6:2FTS    |             | 428.94       | 81.00        | 7.23                 |
| Perfluoro-( <sup>13</sup> C <sub>8</sub> )octanoic acid   | M8PFOA       |             | 420.95       | 376.00       | 7.27                 |
| Perfluoro-(1,2,3,4- <sup>13</sup> C <sub>4</sub> )octanesulfonic acid   | MPFOS        |             | 503.00       | 99.00        | 7.78                 |
| Sodium perfluoro-( <sup>13</sup> C <sub>8</sub> )octanesulfonate  | M8PFOS       |             | 506.91       | 98.90        | 7.89                 |
| Perfluoro-( <sup>13</sup> C <sub>9</sub> )nonanoic acid   | M9PFNA       |             | 471.94       | 427.00       | 7.92                 |
| Perfluoro-(1,2,3,4,5,6- <sup>13</sup> C <sub>6</sub> )decanoic acid   | M6PFDA       |             | 518.92       | 474.00       | 8.49                 |
| Sodium 1 <i>H</i> ,1 <i>H</i> ,2 <i>H</i> ,2 <i>H</i> -perfluoro(1,2- <sup>13</sup> C <sub>2</sub> )decanesulfonate | M2-8:2FTS    |             | 528.94       | 80.90        | 8.50                 |
| Perfluoro-(1,2,3,4,5,6,7- <sup>13</sup> C <sub>7</sub> )undecanoic acid   | M7PFUdA      |             | 569.95       | 525.00       | 8.95                 |

Table 2: MRM transitions and retention times of native PFAS and isotopically labeled PFAS analytical standards.

# Analysis of PFAS in Aqueous Samples with CHROMABOND WAX according to EPA Method 533

## Chromatograms




Figure 1: Chromatogram of a standard solution (concentration,  $\beta = 0.5$  ng/mL)

## Recovery rates




Figure 2: Recovery rate from water sample (concentration,  $\beta = 1$  ng/L,  $n=5$ )

| Analyte                                  | Abbreviation | Recovery rate (%) $\pm$ RSD (%) |
|--|--------------|---------------------------------|
| Perfluoro-3-methoxypropanoic acid        | PFMPA        | 105 $\pm$ 2                     |
| Perfluoro- <i>n</i> -butanoic acid       | PFBA         | 99 $\pm$ 4                      |
| Perfluoro-4-methoxybutanoic acid         | PFMBA        | 95 $\pm$ 3                      |
| Perfluoro- <i>n</i> -pentanoic acid      | PFPeA        | 100 $\pm$ 2                     |
| Perfluoro(2-ethoxyethane) sulfonic acid  | PFEESA       | 94 $\pm$ 3                      |
| Perfluoro- <i>n</i> -butanesulfonic acid | PFBS         | 83 $\pm$ 4                      |

| Analyte   | Abbreviation | Recovery rate (%) $\pm$ RSD (%) |
|---|--------------|---------------------------------|
| Nonafluoro-3,6-dioxaheptanoic acid  | NFDHA        | 84 $\pm$ 4                      |
| 1 <i>H</i> ,1 <i>H</i> , 2 <i>H</i> , 2 <i>H</i> -Perfluorohexane sulfonic acid | 4:2FTS       | 73 $\pm$ 5                      |
| Perfluoro- <i>n</i> -hexanoic acid  | PFHxA        | 95 $\pm$ 3                      |
| Perfluoropentansulfonic acid  | PFPeS        | 69 $\pm$ 6                      |
| Hexafluoropropylene oxide dimer acid  | HFPO-DA      | 100 $\pm$ 5                     |

# Analysis of PFAS in Aqueous Samples with CHROMABOND WAX according to EPA Method 533

| Analyte   | Abbreviation | Recovery rate (%) ± RSD (%) |
|---|--------------|-----------------------------|
| Perfluoro- <i>n</i> -heptanoic acid   | PFHpA        | 95 ± 3                      |
| Perfluoro- <i>n</i> -hexanesulfonic acid  | PFHxS        | 87 ± 4                      |
| 4,8-Dioxa-3 <i>H</i> -perfluorononanoic acid  | ADONA        | 82 ± 2                      |
| 1 <i>H</i> ,1 <i>H</i> , 2 <i>H</i> , 2 <i>H</i> -Perfluorooctane sulfonic acid                                     | 6:2FTS       | 74 ± 9                      |
| Perfluoro- <i>n</i> -heptanesulfonic acid   | PFHpS        | 95 ± 4                      |
| Perfluoro- <i>n</i> -octanoic acid  | PFOA         | 110 ± 10                    |
| Perfluoro- <i>n</i> -octanesulfonic acid  | PFOS         | 83 ± 11                     |
| Perfluoro- <i>n</i> -nonanoic acid  | PFNA         | 94 ± 15                     |
| 9-Chlorohexadecafluoro-3-oxanonane-1-sulfonic acid  | 9Cl-PF3ONS   | 70 ± 9                      |
| Perfluoro- <i>n</i> -decanoic acid  | PFDA         | 80 ± 21                     |
| 1 <i>H</i> ,1 <i>H</i> , 2 <i>H</i> , 2 <i>H</i> -Perfluorodecane sulfonic acid                                     | 8:2FTS       | 75 ± 22                     |
| Perfluoro- <i>n</i> -undecanoic acid  | PFUnDA       | 92 ± 30                     |
| 11-Chloroeicosafluoro-3-oxaundecane-1-sulfonic acid   | 11Cl-PF3OUdS | 61 ± 15                     |
| Perfluoro- <i>n</i> -dodecanoic acid  | PFDoDA       | 72 ± 28                     |
| Surrogates  |              |                             |
| Perfluoro-(2,3,4- <sup>13</sup> C <sub>3</sub> )butanoic acid   | M3PFBA       | 86 ± 2                      |
| Perfluoro-( <sup>13</sup> C <sub>4</sub> )butanoic acid   | M4PFBA       | 85 ± 4                      |
| Sodium perfluoro-(2,3,4- <sup>13</sup> C <sub>3</sub> )butanesulfonate  | M3PFBS       | 86 ± 3                      |
| Sodium 1 <i>H</i> ,1 <i>H</i> ,2 <i>H</i> ,2 <i>H</i> -perfluoro(1,2- <sup>13</sup> C <sub>2</sub> )hexanesulfonate | M2-4:2FTS    | 84 ± 9                      |
| Tetrafluoro-2-heptafluoropropoxy- <sup>13</sup> C <sub>3</sub> -propanoic acid                                      | M3HPFO-DA    | 99 ± 3                      |
| Perfluoro-(1,2,3,4- <sup>13</sup> C <sub>4</sub> )heptanoic acid  | M4PFHpA      | 100 ± 2                     |
| Sodium perfluoro-(1,2,3- <sup>13</sup> C <sub>3</sub> )hexanesulfonate  | M3PFHxS      | 84 ± 5                      |
| Perfluoro-(1,2,3,4- <sup>13</sup> C <sub>4</sub> )octanoic acid   | MPFOA        | 111 ± 5                     |
| Sodium 1 <i>H</i> ,1 <i>H</i> ,2 <i>H</i> ,2 <i>H</i> -perfluoro(1,2- <sup>13</sup> C <sub>2</sub> )octanesulfonate | M2-6:2FTS    | 127 ± 25                    |
| Perfluoro-( <sup>13</sup> C <sub>8</sub> )octanoic acid   | M8PFOA       | 97 ± 3                      |

| Analyte  | Abbreviation | Recovery rate (%) ± RSD (%) |
|--|--------------|-----------------------------|
| Perfluoro-(1,2,3,4- <sup>13</sup> C <sub>4</sub> )octanesulfonic acid  | MPFOS        | 80 ± 3                      |
| Sodium perfluoro-( <sup>13</sup> C <sub>8</sub> )octanesulfonate   | M8PFOS       | 92 ± 2                      |
| Perfluoro-( <sup>13</sup> C <sub>9</sub> )nonanoic acid  | M9PFNA       | 89 ± 1                      |
| Perfluoro-(1,2,3,4,5,6- <sup>13</sup> C <sub>6</sub> )decanoic acid  | M6PFDA       | 84 ± 2                      |
| Sodium 1 <i>H</i> ,1 <i>H</i> ,2 <i>H</i> ,2 <i>H</i> -perfluoro(1,2- <sup>13</sup> C <sub>2</sub> )decansulfonate | M2-8:2FTS    | 94 ± 20                     |
| Perfluoro-(1,2,3,4,5,6,7- <sup>13</sup> C <sub>7</sub> )undecanoic acid  | M7PFUdA      | 80 ± 5                      |

Table 3: Recovery rates for the presented SPE method using CHROMABOND® WAX, 150 mg, 6 mL (n=5).


## Conclusion

This application note presents the reliable and successful determination of 25 PFAS according to EPA method 533 from drinking water. By using the SPE column, CHROMABOND® WAX, it was possible to achieve high recovery rates with good reproducibility. CHROMABOND® WAX was optimized for PFAS analysis and provides various interaction types like ionic, hydrophobic, hydrogen bonds and dipole-dipole interactions for the enrichment of a broad spectrum of PFAS. The sorbent is specially recommended for PFAS analysis because of its very low blind value levels. Most of the PFAS show recovery rates between 80 % to 110 %. The identification and the quantification of PFAS in food were finally carried out by ESI mass spectrometry on a NUCLEODUR® PFAS column.

## References

[1] METHOD 533: DETERMINATION OF PER- AND POLYFLUOROALKYL SUBSTANCES IN DRINKING WATER BY ISOTOPE DILUTION ANION EXCHANGE SOLID PHASE EXTRACTION AND LIQUID CHROMATOGRAPHY/TANDEM MASS SPECTROMETRY; December 2019.

MACHEREY NAGEL shall not be liable for errors contained herein or for incidental or consequential damages in connection with the furnishing, performance, or use of this material. Information, descriptions, and specifications in this publication are subject to change without notice.



### AUTHORIZED DISTRIBUTOR

MZ-Analysentechnik GmbH, Barcelona-Allee 17 • D-55129 Mainz  
Tel +49 6131 880 96-0, Fax +49 6131 880 96-20  
e-mail: info@mz-at.de, www.mz-at.de

www.mn-net.com

# MACHEREY-NAGEL



MACHEREY-NAGEL GmbH & Co. KG  
Valenciennner Str. 11  
52355 Düren · Germany

DE Tel.: +49 24 21 969-0 info@mn-net.com  
CH Tel.: +41 62 388 55 00 sales-ch@mn-net.com  
FR Tel.: +33 388 68 22 68 sales-fr@mn-net.com  
US Tel.: +1 888 321 62 24 sales-us@mn-net.com

