

Chromatography

# CHROMAFIL® Xtra



Pure Filters

**MACHERY-NAGEL**

[www.mn-net.com](http://www.mn-net.com)



*Since 1911*



# CHROMAFIL® - Disposable filters from MACHEREY-NAGEL



## Contents

➤ Introduction	3
➤ Technical information	5-9
➤ Specifications	11
➤ Quality control	12
➤ Packages and membranes	13-14
➤ Chemical compatibility	15
➤ Optimal use of CHROMAFIL®	16
➤ How to select the optimal CHROMAFIL®	17
➤ Ordering information	18-29
➤ Technical support and customer service	31



## CHROMAFIL®

Ideal for GC, HPLC and UHPLC sample clarification

### Introduction

CHROMAFIL® syringe filters are used for filtration of suspended matter from liquid samples (1–100 mL) or gases. With CHROMAFIL®, rapid purification and removal of particles is very simple: just place the filter on the syringe, and you are ready for filtration. Special manipulations are not required. Contamination of sensitive instrumentation by solid impurities can be avoided, thus increasing lifetime of chromatographic columns and equipment. The filter can be used for the sample preparation for HPLC, GC, ICP, AAS, TOC, DOC, IR, NMR, photometry, spectroscopy, . . .

- different membrane types to meet multiple filtration needs
- low content of extractable compounds ensure reliable analyses
- superior chromatography column protection helps extend column life
- fast flow geometry for easy filtration
- low hold-up volume for maximum filtrate recovery
- HPLC certified
- designed to be compatible for use on all common automated filtration systems, e.g. SOTAX dissolution systems


### CHROMAFIL® Xtra

labeled for method validation and certification

**Xtra** imprint for direct identification of the membrane type, diameter and pore size

**Xtra** low bleeding PP housing

**Xtra** color-free plain polypropylene







## Technical Information

◆ **Low content of extractable substances** due to a high density polypropylene housing combined with ultrapure filtration membranes.

In comparison to filters made of polycarbonate, polyacrylate or polystyrene, all CHROMAFIL® filters are resistant against nearly all organic solvents.

(see list of chemical compatibility on page 15)

## HPLC-test

Conditions: 2 mL of the solvent (specified on top of the chromatograms) were applied to the filter; 100 µL of the filtrate were injected to the HPLC.

Eluent A: water

Eluent B: acetonitrile

Gradient: 10 % → 95 % B in 25 minutes


Flow rate: 0.5 mL/min

Sensitivity: -5 to 10 mAU at 254 nm


Column: 125/4 NUCLEODUR® C<sub>18</sub> Gravity 5 µm  
(REF: 760100.40)

### Acetonitrile:

Eluate:



Blanc:







## CHROMAFIL® - Disposable filters from MACHERY-NAGEL

### Low Bleeding PP housing

Even a treatment with very aggressive solvents/reagents does not lead to significant blind values or extractables.

For a proof, a filter was rinsed with 0.5 mL MSTFA (N-methyl-N-trimethylsilyl-trifluoroacetamide), a very powerful silylation reagent. The result is shown in the GC/MS chromatogram.





To provide the **lowest content of extractable substances**, the housing of every CHROMAFIL® filter is ultrasonically sealed

- The filters are **welded, not glued**, because glue may have extractable ingredients
- The welding leads to a tight connection between the both parts, thus the filter can be used in both directions. No fluid can leak from the filter housing.

The special **thick rim** of the housing is ideal for use in laboratory robots (e. g. SOTAX, Benchmate™, . . .).



For a **safe connection on the „high pressure“ side** every CHROMAFIL® filter provides a **Luer lock on inlet**.



## Luer outlet

- For the 3, 13 and 25 mm diameter filter: standard luer outlet  
This luer configuration offers low hold-up volume and easy filtration into autosampler vials and NMR tubes

Filter inlet and filter outlet can be fitted to all CHROMABOND® columns and accessories for selective sample preparation with the aid of a special adaptor.



## CHROMAFIL® - Disposable filters from MACHERY-NAGEL



### No breakage of the membrane due to a stabilizing “crash” plate

- The sample fluid is deviated in four lanes by the “crash” plate and does not directly hit the membrane. The resulting pressure distribution protects the membrane against breakage.

### Optimal flow geometry by star-shaped distribution plate


- The fluid penetrates the membrane on the whole surface, not only on a small area; the **filter will not clog rapidly**, which guarantees in a high flow efficiency.



### Different pore sizes for multiple filtration application

Available **pore sizes** 0.2 and 0.45  $\mu\text{m}$  (additional: PET filters with 1.2  $\mu\text{m}$ , glass fiber filters with 1  $\mu\text{m}$ , PES filters with 5  $\mu\text{m}$ ). Filters with 0.45  $\mu\text{m}$  pore size remove fine particles which can clog chromatography columns. **0.2  $\mu\text{m}$  pore size filters are recommended for filtration of UHPLC samples.**






## Filter Sizes

3, 13 and 25 mm effective membrane diameter. The small diameter filters are especially recommended for very small samples, which require extremely low dead volumes.


Sample volume	Recommended membrane diameter	Dead volume	Filtration area
≤ 1 mL	3 mm	5 µL	0.07 cm <sup>2</sup>
1-5 mL	13 mm	30 µL	1.33 cm <sup>2</sup>
5-100 mL	25 mm	80 µL	4.91 cm <sup>2</sup>



All filters can be **autoclaved** at **121 °C** and **1.1 bar** for 30 min.

25 mm CHROMAFIL® filters are designed to be **100% compatible and reliable for use with the SOTAX AT70 smart** fully automated dissolution testing systems.





## A specification data sheet

is available for all membranes and filter diameters

## Enhanced quality control for better results

MN certifies CHROMAFIL® syringe filters to be low in UV absorbing extractables.

All filters and membrane types have been HPLC tested for compatibility with the most common HPLC solvents (methanol, water, acetonitrile, see test chromatograms).

HPLC-test certificates are available for every membrane type.

Please visit: [www.mn-net.com](http://www.mn-net.com)

Example of specification data sheet



## CHROMAFIL® - Disposable filters from MACHEREY-NAGEL



### Pressure stability of CHROMAFIL® syringe filter housing 12 bar

The „blue“ test: membrane, pressure and filtration batch test with blue colored silica particles in matching particle sizes provides an excellent method to find leaks or membrane deviations.



## Package sizes

packs of 100 or 400 (BigBoxes) for 25 mm Ø filter  
packs of 100 for 13 mm Ø filter  
packs of 100 for 3 mm Ø filter  
packs of 50 for sterile filter

## Different membrane materials for multiple filtration applications


Depending on your filtration task you can choose filter membranes made from different materials:

<b>Polyester</b>	(PET) with or without <b>glass fiber prefilter*</b>
<b>Regenerated cellulose</b>	(RC) with or without <b>glass fiber prefilter*</b>
<b>Polytetrafluoroethylene</b>	(PTFE)
<b>Hydrophilized Polytetrafluoroethylene</b>	(H-PTFE)
<b>Cellulose mixed esters</b>	(MV)
<b>Cellulose acetate</b>	(CA) sterile and non-sterile
<b>Polyamide / Nylon</b>	(PA)
<b>Polyethersulfone</b>	(PES)
<b>Polyvinylidene difluoride</b>	(PVDF) with or without <b>glass fiber prefilter*</b>
<b>Glass fiber</b>	(GF)

\* Filters with (nom. 1 µm) GF prefilter provide a 2-4 times greater throughput than filter without prefilter for extremely viscous and most difficult-to-filter samples



## CHROMAFIL® - Disposable filters from MACHEREY-NAGEL



### Combi syringe filters with a coarse glass fiber (GF) prefilter and a small-pore membrane as main filter

#### User benefits:

- for solutions with a high load of particulate matter: lower back pressure, easy filtration
- for high yields of filtrate: more mL of pure filtrate per filter

#### The technology:

The glass fiber membrane (1  $\mu\text{m}$ ) removes coarse particles, before they can block the fine main membrane. This results in a better filtration efficiency, particularly for highly contaminated samples.

Housing: solvent-resistant, ultra low bleed polypropylene

Inlet: Luer lock

Outlet: Luer

Pore diameter: 1.0/0.20  $\mu\text{m}$  or 1.0/0.45  $\mu\text{m}$

Filter diameter: 25 mm

Void volume: < 80  $\mu\text{L}$

Packing unit: 100 filters / BigBoxes with 400 filters

#### Available membranes with GF-prefilter:

Polyester (PET)

Regenerated cellulose (RC)

Polyvinylidene Difluoride (PVDF)

# CHROMAFIL® - Disposable filters from MACHERY-NAGEL



Solvent	Material										
	MV	CA	RC	PA	PTFE	H-PTFE	PVDF	PES	PET	GF	PP
Acetaldehyde	⊖	⊖	⊕	⊕	⊕	⊕	⊕	⊕	⊕	⊕	⊕
Acetic acid, 100 %	⊖	⊖	⊖	⊖	⊕	⊕	⊕	⊕	⊕	⊕	⊕
Acetone	⊖	⊖	⊕	⊕	⊕	⊕	⊖	⊖	⊕	⊕	⊕
Acetonitrile	⊖	⊖	⊕	⊕	⊕	⊕	⊕	⊕	⊕	⊕	⊕
Ammonia, 25 %	⊖	⊖	⊕	⊖	⊕	⊕	⊕	⊕	⊕	⊕	⊕
Benzene	⊕	⊕	⊕	⊕	⊕	⊕	⊕	⊕	⊕	⊕	⊕
n-Butanol	⊕	⊕	⊕	⊕	⊕	⊕	⊕	⊕	⊕	⊕	⊕
Cyclohexane	⊕	⊕	⊕	⊕	⊕	⊕	⊕	⊕	⊕	⊕	⊕
Dichloromethane	⊕	⊕	⊖	⊖	⊕	⊕	⊕	⊖	⊕	⊕	⊖
Diethyl ether	⊕	⊕	⊕	⊕	⊕	⊕	⊕	⊕	⊕	⊕	⊕
Dimethylformamide	⊖	⊖	⊕	⊕	⊕	⊕	⊖	⊖	⊕	⊕	⊕
1,4-Dioxane	⊖	⊖	⊕	⊕	⊕	⊕	⊕	⊖	⊕	⊕	⊕
Ethanol	⊖	⊕	⊕	⊕	⊕	⊕	⊕	⊕	⊕	⊕	⊕
Ethyl acetate	⊖	⊖	⊕	⊕	⊕	⊕	⊕	⊕	⊕	⊕	⊕
Ethylene glycol	⊕	⊕	⊕	⊕	⊕	⊕	⊕	⊕	⊕	⊕	⊕
Formic acid, 100 %	⊕	⊖	⊕	⊖	⊕	⊕	⊕	⊕	⊕	⊕	⊕
Hydrochloric acid, 30 %	⊖	⊖	⊖	⊖	⊕	⊕	⊕	⊕	⊕	⊕	⊕
Methanol	⊖	⊖	⊕	⊕	⊕	⊕	⊕	⊕	⊕	⊕	⊕
Nitric acid, 65 %	⊖	⊖	⊖	⊖	⊕	⊕	⊕	⊕	⊕	⊕	⊖
Oxalic acid, 10 % aqueous	⊕	⊖	⊕	⊖	⊕	⊕	⊕	⊕	⊕	⊕	⊕
Petroleum ether	⊕	⊕	⊕	⊕	⊕	⊕	⊕	⊕	⊕	⊕	⊕
Phosphoric acid, 80 %	⊖	⊖	⊕	⊖	⊕	⊕	⊕	⊕	⊕	⊕	⊕
Potassium hydroxide, 1 mol/L	⊖	⊖	⊕	⊕	⊕	⊕	⊕	⊕	⊕	⊕	⊕
2-Propanol	⊕	⊕	⊕	⊕	⊕	⊕	⊕	⊕	⊕	⊕	⊕
Sodium hydroxide, 1 mol/L	⊖	⊖	⊕	⊕	⊕	⊕	⊕	⊕	⊕	⊕	⊕
Tetrachloromethane	⊕	⊖	⊕	⊕	⊕	⊕	⊕	⊕	⊕	⊕	⊕
Tetrahydrofuran	⊖	⊖	⊕	⊕	⊕	⊕	⊕	⊖	⊕	⊕	⊕
Toluene	⊕	⊖	⊕	⊕	⊕	⊕	⊕	⊕	⊕	⊕	⊕
Trichloroethylene	⊕	⊕	⊕	⊕	⊕	⊕	⊕	⊕	⊕	⊕	⊕
Trichloromethane (Chloroform)	⊕	⊖	⊕	⊖	⊕	⊕	⊕	⊖	⊕	⊕	⊖
Urea	⊕	⊕	⊕	⊕	⊕	⊕	⊕	⊕	⊕	⊕	⊕
Water	⊕	⊕	⊕	⊕	⊕	⊕	⊕	⊕	⊕	⊕	⊕
Xylene	⊕	⊕	⊕	⊕	⊕	⊕	⊕	⊕	⊕	⊕	⊕

The table lists the chemical compatibility of our CHROMAFIL® materials. The chemical compatibility depends on several parameters such as time, pressure, temperature and concentration. In most cases, CHROMAFIL® filters will have only short contact with a solvent. In these cases they may be used despite of limited compatibility. For example, a PTFE filter with PP housing does not release any UV-detectable substances during filtration of 5 mL THF, although PP shows only limited resistance towards THF.

## Data not guaranteed.

⊕ resistant, ⊖ not resistant, ⊕ limited resistance

MV = cellulose mixed esters

CA = cellulose acetate

RC = regenerated cellulose

PA = polyamide (Nylon)

PTFE = polytetrafluoroethylene

H-PTFE = hydrophilized polytetrafluoroethylene

PVDF = polyvinylidene difluoride

PES = polyethersulfone

PET = polyester

GF = glass fiber

PP = polypropylene (housing material)



# CHROMAFIL® - Disposable filters from MACHEREY-NAGEL

## Optimal use of CHROMAFIL® syringe filter

For achieving the full benefits of filtration we recommend the following instructions.



Draw up the sample into the syringe. Then draw approximately 1 mL of air into the syringe. The air helps to minimize the remaining fluid in the filter.

Plug the CHROMAFIL® syringe filter onto the syringe with the luer connection. Ensure a tight connection by gently turning.

Start with gentle pressure to filter your sample into a vial\*. This helps to assure maximum throughput.

### Tips / additional information

We recommend either discarding the first 1 mL or rinsing the filter unit with 1 mL of primary solvent before sample filtration.

In order to avoid breakage of the membrane only syringes with volumes of 10 mL or higher should be used.

Do not reuse syringe filters  
Do not use at temperatures above 55 °C (131 °F)











**Warning:** CHROMAFIL® syringe filters are intended for laboratory use only. Do not use CHROMAFIL® syringe filters for direct patient care applications.

\* MACHEREY-NAGEL offers a wide range of vials and caps. More information at [www.mn-net.com/vials](http://www.mn-net.com/vials)





## How to select the optimal CHROMAFIL® syringe filter

1. Filter size		2. Pore size of filter membrane	3. Membrane type
<p>sample volume</p>  <p>5-100 mL</p>	<p>filter size</p>  <p>25 mm</p>	<p>pore size</p>  <p>0.45 µm</p>	<p>properties of sample</p> <p><b>aqueous, polar hydrophilic</b> low particle-load</p> <p><b>high particle-load, prefiltration required</b></p> <p><b>mid-polar e.g. HPLC eluents</b></p> <p><b>proteins</b> protein needed</p> <p><b>remove protein</b></p> <p><b>strong acids and bases</b></p> <p><b>organic, nonpolar, hydrophobic</b> low particle-load</p> <p><b>high particle-load, prefiltration required</b></p>
<p>1-5 mL</p> 	 <p>13 mm</p>	 <p>0.20 µm</p>	<p><b>PET</b></p> <p><b>H-PTFE</b></p> <p><b>MV</b></p> <p><b>RC</b></p> <p><b>GF/PET</b></p> <p><b>GF/RC</b></p> <p><b>GF/PVDF</b></p> <p><b>PET</b></p> <p><b>PA</b></p> <p><b>RC</b></p> <p><b>CA</b></p> <p><b>PVDF</b></p> <p><b>PES</b></p> <p><b>GF</b></p> <p><b>GF/PET</b></p> <p><b>GF/PVDF</b></p> <p><b>H-PTFE</b></p> <p><b>PTFE</b></p> <p><b>PTFE</b></p> <p><b>PET</b></p> <p><b>GF/PET</b></p> <p><b>GF/PVDF</b></p>
<p>&lt; 1 mL</p> 	 <p>3 mm</p>		



## CHROMAFIL® - Disposable filters from MACHEREY-NAGEL



### Polyester (PET)

- ◆ hydrophilic multipurpose membrane
- ◆ for polar as well as nonpolar solvents  
**the HPLC filter**, especially suited for mixtures of water and organic solvents for TOC/DOC determination, not cytotoxic, does not inhibit the growth of microorganisms and higher cells
- ◆ PET filters with integrated glass fiber prefilter (**GF/PET**) are recommended for solutions with a high load of particulate matter or for highly viscous solutions

### Ordering information · CHROMAFIL®




	Type	Pore size [µm]	Membrane diameter [mm]	Color code		Standard pack		BIG-BOX	
				top	bottom	filters / pack	REF	filters / pack	REF
	Xtra PET-20/25	0.20	25	labeled		100	<b>729221</b>	400	<b>729221.400</b>
	Xtra PET-45/25	0.45	25	labeled		100	<b>729220</b>	400	<b>729220.400</b>
	Xtra PET-120/25	1.2	25	labeled		100	<b>729229</b>	400	<b>729229.400</b>
	Xtra PET-20/13	0.20	13	labeled		100	<b>729222</b>		
	Xtra PET-45/13	0.45	13	labeled		100	<b>729223</b>		
<b>Combi Filters</b>									
	GF/PET-20/25	1.0/0.20	25	blue	orange	100	<b>729032</b>	400	<b>729032.400</b>
	GF/PET-45/25	1.0/0.45	25	black	orange	100	<b>729033</b>	400	<b>729033.400</b>



## Regenerated Cellulose (RC)

- ◆ hydrophilic membrane with very low adsorption
- ◆ for aqueous and organic / aqueous liquids, i. e. polar and medium polar sample solutions
- ◆ binding capacity for proteins 84 µg per 25 mm filter
- ◆ RC filters with integrated glass fiber prefilter\* (**GF / RC**) are recommended for solutions with a high load of particulate matter or for highly viscous solutions

## Ordering information · CHROMAFIL®

	Type	Pore size [µm]	Membrane diameter [mm]	Color code		Standard pack		BIG-BOX	
				top	bottom	filters / pack	REF	filters/pack	REF
	Xtra RC-20/25	0.20	25	labeled		100	<b>729230</b>	400	<b>729230.400</b>
	Xtra RC-45/25	0.45	25	labeled		100	<b>729231</b>	400	<b>729231.400</b>
	Xtra RC-20/13	0.20	13	labeled		100	<b>729236</b>		
	Xtra RC-45/13	0.45	13	labeled		100	<b>729237</b>		
<b>Combi filters</b>									
	GF/RC-20/25	1.0/0.20	25	blue	blue	100	<b>729050</b>	400	<b>729050.400</b>
	GF/RC-45/25	1.0/0.45	25	black	blue	100	<b>729051</b>	400	<b>729051.400</b>

\* glass fiber exhibits a high protein-binding capacity



## CHROMAFIL® - Disposable filters from MACHERY-NAGEL



### Polytetrafluoroethylene (PTFE)

- ◆ hydrophobic membrane
- ◆ for nonpolar liquids and gases
- ◆ very resistant towards all kinds of solvents as well as acids and bases  
flushing with alcohol, followed by water, makes the originally hydrophobic membrane more hydrophilic

### Ordering information · CHROMAFIL®


	Type	Pore size [µm]	Membrane diameter [mm]		Standard pack		BIG-BOX	
					filters / pack	REF	filters / pack	REF
	Xtra PTFE-20/25	0.20	25	labeled	100	<b>729207</b>	400	<b>729207.400</b>
	Xtra PTFE-45/25	0.45	25	labeled	100	<b>729205</b>	400	<b>729205.400</b>
	Xtra PTFE-20/13	0.20	13	labeled	100	<b>729208</b>		
	Xtra PTFE-45/13	0.45	13	labeled	100	<b>729209</b>		
	O-20/3	0.20	3		100	<b>729014</b>		
	O-45/3	0.45	3		100	<b>729015</b>		



## Hydrophilized polytetrafluoroethylene (H-PTFE)

- ◆ hydrophobic membrane with additional hydrophilic properties
- ◆ for polar and nonpolar sample solutions
- ◆ resistant towards all kinds of solvents as well as acids and bases

## Ordering information · CHROMAFIL®

	Type	Pore size [µm]	Membrane diameter [mm]		Standard pack		BIG-BOX	
					filters / pack	REF	filters / pack	REF
	Xtra H-PTFE-20/25	0.20	25	labeled	100	<b>729245</b>	400	<b>729245.400</b>
	Xtra H-PTFE-45/25	0.45	25	labeled	100	<b>729246</b>	400	<b>729246.400</b>
	Xtra H-PTFE-20/13	0.20	13	labeled	100	<b>729256</b>		
	Xtra H-PTFE-20/13	0.45	13	labeled	100	<b>729257</b>		



## CHROMAFIL® - Disposable filters from MACHEREY-NAGEL



### Cellulose Mixed Ester (MV)

- ◆ hydrophilic membrane with very low adsorption
- ◆ for aqueous or polar solutions

### Ordering information · CHROMAFIL®

	Type	Pore size [µm]	Membrane diameter [mm]		Standard pack		BIG-BOX	
					filters / pack	REF	filters / pack	REF
	Xtra MV-20/25	0.20	25	labeled	100	<b>729206</b>	400	<b>729206.400</b>
	Xtra MV-45/25	0.45	25	labeled	100	<b>729204</b>	400	<b>729204.400</b>



## Cellulose Acetate (CA)

- ◆ hydrophilic membrane
- ◆ for filtration of water-soluble oligomers and polymers, especially suited for biological macromolecules
- ◆ very high stability in aqueous solutions
- ◆ binding capacity for proteins 21 µg per 25 mm filter
- ◆ also available in a sterile package (S) for filtration under sterile conditions (each filter individually sealed)

## Ordering information · CHROMAFIL®

	Type	Pore size [µm]	Membrane diameter [mm]	Color code		Standard pack		BIG-BOX	
				top	bottom	filters / pack	REF	filters / pack	REF
	Xtra CA-20/25	0.20	25	labeled		100	<b>729226</b>	400	<b>729226.400</b>
	Xtra CA-45/25	0.45	25	labeled		100	<b>729227</b>	400	<b>729227.400</b>
	Xtra CA-20/13	0.20	13	labeled		100	<b>729254</b>		
	Xtra CA-45/13	0.45	13	labeled		100	<b>729255</b>		
<b>Sterile filters</b>									
	CA-20/25 (S)	0.20	25	yellow	red	50	<b>729024</b>		
	CA-45/25 (S)	0.45	25	colorless	red	50	<b>729025</b>		



## CHROMAFIL® - Disposable filters from MACHERY-NAGEL



### Polyamide (PA) = Nylon

- moderately hydrophilic membrane
- for aqueous and organic / aqueous medium polar liquids

### Ordering information · CHROMAFIL®

	Type	Pore size [µm]	Membrane diameter [mm]		Standard pack		BIG-BOX	
					filters / pack	REF	filters / pack	REF
	Xtra PA-20/25	0.20	25	labeled	100	<b>729212</b>	400	<b>729212.400</b>
	Xtra PA-45/25	0.45	25	labeled	100	<b>729213</b>	400	<b>729213.400</b>
	Xtra PA-20/13	0.20	13	labeled	100	<b>729248</b>		
	Xtra PA-45/13	0.45	13	labeled	100	<b>729249</b>		
	AO-20/3	0.20	3		100	<b>729010</b>		
	AO-45/3	0.45	3		100	<b>729011</b>		





## Polyethersulfone (PES)

- ◆ hydrophilic membrane
- ◆ for aqueous and slightly organic liquids with higher flow rates
- ◆ very low adsorption for pharmaceuticals and proteins
- ◆ good stability against organic acids and bases
- ◆ binding capacity for proteins 29 µg per 25 mm filter

## Ordering information · CHROMAFIL®

	Type	Pore size [µm]	Membrane diameter [mm]		Standard pack		BIG-BOX	
					filters / pack	REF	filters / pack	REF
	Xtra PES-20/25	0.20	25	labeled	100	<b>729240</b>	400	<b>729240.400</b>
	Xtra PES-45/25	0.45	25	labeled	100	<b>729241</b>	400	<b>729241.400</b>
	Xtra PES-500/25	5.0	25	labeled	100	<b>729242</b>	400	<b>729242.400</b>



## CHROMAFIL® - Disposable filters from MACHEREY-NAGEL



### Polyvinylidene Difluoride (PVDF)

- hydrophilic membrane
- for aqueous solutions, water-soluble oligomers and polymers like proteins
- low binding capacity for proteins 20 µg per 25 mm filter
- PVDF filters with integrated glass fiber prefilter\* (**GF / PVDF**) are recommended for filtration of biological samples with high particle loads.

### Ordering information · CHROMAFIL®

	Type	Pore size [µm]	Membrane diameter [mm]	Color code		Standard pack		BIG-BOX	
				top	bottom	filters / pack	REF	filters / pack	REF
	Xtra PVDF-20/25	0.20	25	labeled		100	<b>729218</b>	400	<b>729218.400</b>
	Xtra PVDF-45/25	0.45	25	labeled		100	<b>729219</b>	400	<b>729219.400</b>
	Xtra PVDF-20/13	0.20	13	labeled		100	<b>729243</b>		
	Xtra PVDF-45/13	0.45	13	labeled		100	<b>729244</b>		
<b>Combi filters</b>									
	GF/PVDF-45/25	1.0/0.45	25	black	white	100	<b>729039</b>	400	<b>729039.400</b>

\* glass fiber exhibits a high protein-binding capacity



## Glass Fiber (GF)

- ◆ inert filter, nominal pore size 1 µm, allows higher flow rates than small pore filters
- ◆ for solutions with high loads of particulate matter or for highly viscous solutions (e. g. soil samples, fermentation broths)
- ◆ as prefilters for other CHROMAFIL® filters, they prevent clogging of the membrane

## Ordering information · CHROMAFIL®

	Type	Pore size [µm]	Membrane diameter [mm]		Standard pack		BIG-BOX	
					filters / pack	REF	filters / pack	REF
	Xtra GF-100/25	nom. 1.0	25	labeled	100	<b>729228</b>	400	<b>729228.400</b>
	GF-100/13	nom. 1.0	13	labeled	100	<b>729234</b>		



## CHROMAFIL® - Disposable filters from MACHEREY-NAGEL




### Filtration cartridges

- Filtration cartridges for sample clarification under vacuum (e.g., using the CHROMABOND® vacuum manifold or SPE automation systems like Gilson Aspec™, Rapidtrace) or by gravity flow.
- Cartridge sizes 3 mL and 6 mL
- Different membranes (PET, RC, PTFE, PVDF, GF) and pore sizes (0.2, 0.45 and 1.0 µm). The membrane materials correspond to the respective CHROMAFIL® syringe filters.

### Ordering information · CHROMAFIL® filtration cartridges

	Type	Pore size [µm]	Pack of	Column volume REF	
				3 mL	6 mL
	PET (polyester)	0.20	100	<b>730578.320</b>	<b>730578.620</b>
	PET (polyester)	0.45	100	<b>730578.345</b>	<b>730578.645</b>
	RC (regenerated cellulose)	0.20	100	<b>730068.320</b>	<b>730068.620</b>
	RC (regenerated cellulose)	0.45	100	<b>730068.345</b>	<b>730068.645</b>
	PTFE (polytetrafluoroethylene)	0.20	100	<b>730570.320</b>	<b>730570.620</b>
	PTFE (polytetrafluoroethylene)	0.45	100	<b>730570.345</b>	<b>730570.645</b>
	PVDF (polyvinylidene difluoride)	0.20	100	<b>730579.320</b>	<b>730579.620</b>
	PVDF (polyvinylidene difluoride)	0.45	100	<b>730579.345</b>	<b>730579.645</b>
GF (glass fiber)	nom. 1.0	100	<b>730517.3100</b>	<b>730517.6100</b>	



## MULTI 96 filter plates

- ◆ 96-well polypropylene plates for simultaneous filtration of 96 samples
- ◆ Advantages of this high-throughput system:
  - Economical by saving time and solvent
  - Use of multi-channel pipettors facilitates liquid transfer steps
  - Readily adaptable to all common automated/robotic handling systems
  - Minimized dead volume ( $\leq 40 \mu\text{L}$ )
- ◆ Membrane materials correspond to the respective CHROMAFIL® syringe filters

## Ordering information · CHROMAFIL® MULTI 96 filter plates

	Description	Pack of	REF
	Filter plates with cellulose mixed ester filter elements (0.20 $\mu\text{m}$ )	1	<b>738770.M</b>
	Filter plates with cellulose mixed ester filter elements (0.45 $\mu\text{m}$ )	1	<b>738771.M</b>
	Filter plates with RC filter elements (regenerated cellulose, 0.20 $\mu\text{m}$ )	1	<b>738656.M</b>
	Filter plates with RC filter elements (regenerated cellulose, 0.45 $\mu\text{m}$ )	1	<b>738657.M</b>
	Filter plates with PTFE filter elements (0.20 $\mu\text{m}$ )	1	<b>738660.M</b>
	Filter plates with PTFE filter elements (0.45 $\mu\text{m}$ )	1	<b>738661.M</b>
	Filter plates with PTFE filter elements (1.0 $\mu\text{m}$ )	1	<b>738662.M</b>
	Filter plates with PTFE filter elements (3.0 $\mu\text{m}$ )	1	<b>738663.M</b>
	Filter plates with PE filter elements (20 $\mu\text{m}$ )	1	<b>738655.M</b>
	Filter plates with PE filter elements (50 $\mu\text{m}$ )	1	<b>738659.M</b>
	Filter plates with glass fiber filter elements (nominal 1 $\mu\text{m}$ )	1	<b>738655.2M</b>
	Filter plates with glass fiber filter elements (nominal 3 $\mu\text{m}$ )	1	<b>738658.M</b>





## Technical Support and Customer Service

### ... we Meet your Needs

If you have any questions concerning CHROMAFIL® filters or our chromatography program, or if you are looking for solutions to a special application, please feel free to contact us:

[www.mn-net.com](http://www.mn-net.com) · [info@mn-net.com](mailto:info@mn-net.com)

Our website offers **more than 3000 applications** which might already solve your analytical questions.

Please visit: [www.mn-net.com/apps](http://www.mn-net.com/apps)



**Germany  
and international**

phone +49 24 21 969-0  
toll-free 0800 2616 000  
fax +49 24 21 969-199 or -198  
e-mail [sales-de@mn-net.com](mailto:sales-de@mn-net.com)



**USA**

phone +1 484 821 0984  
toll-free 888-321-6224 (MACH)  
fax +1 484 821 1272  
e-mail [sales-us@mn-net.com](mailto:sales-us@mn-net.com)



**France**

phone +33 388 68 22 68  
fax +33 388 51 76 88  
e-mail [sales-fr@mn-net.com](mailto:sales-fr@mn-net.com)



**Switzerland**

phone +41 62 388 55 00  
fax +41 62 388 55 05  
e-mail [sales-ch@mn-net.com](mailto:sales-ch@mn-net.com)

**“FilterFinder” online:  
THE cross reference for syringe filter  
[www.mn-net.com/filterfinder](http://www.mn-net.com/filterfinder)**





HPLC



GC



TLC




SPE and Flash




Syringe filters



Vials and caps



... we Meet your Needs



**AUTHORIZED DISTRIBUTOR**

MZ-Analysentechnik GmbH, Barcelona-Allee 17• D-55129 Mainz

Tel +49 6131 880 96-0, Fax +49 6131 880 96-20

e-mail: info@mz-at.de, www.mz-at.de

[www.mn-net.com](http://www.mn-net.com)

**MACHERY-NAGEL**



**MACHERY-NAGEL GmbH & Co. KG** · Neumann-Neander-Str. 6–8 · 52355 Düren · Germany

Germany

and international:

Tel.: +49 24 21 969-0

Fax: +49 24 21 969-199

E-mail: info@mn-net.com

Switzerland:

**MACHERY-NAGEL AG**

Tel.: +41 62 388 55 00

Fax: +41 62 388 55 05

E-mail: sales-ch@mn-net.com

France:

**MACHERY-NAGEL EURL**

Tel.: +33 388 68 22 68

Fax: +33 388 51 76 88

E-mail: sales-fr@mn-net.com

USA:

**MACHERY-NAGEL Inc.**

Tel.: +1 484 821 0984

Fax: +1 484 821 1272

E-mail: sales-us@mn-net.com



Since 1911