

# How to Choose a Column

# Contents

<b>Particle Technology</b> .....	<b><a href="#">85</a></b>
Ethylene Bridged Hybrid (BEH) Particle Technology.....	<a href="#">85</a>
Charged Surface Hybrid (CSH) Particle Technology.....	<a href="#">86</a>
High Strength Silica (HSS) Particle Technology.....	<a href="#">86</a>
Solid-Core Particle Technology.....	<a href="#">86</a>
Reversed Phase and HILIC Chemistries.....	<a href="#">86</a>
<b>Column Selection</b> .....	<b><a href="#">87</a></b>
USP "L" Column Listing.....	<a href="#">87</a>
<b>Column Configurations for Any LC System</b> .....	<b><a href="#">91</a></b>
Column Nomenclature.....	<a href="#">91</a>
Column Configuration.....	<a href="#">91</a>
L/dp Comparison Chart for LC Columns.....	<a href="#">92</a>

# How to Choose a Column

Separation scientists continue to search for innovative solutions to improve chromatographic performance. With a wide array of column choices and formats, they have the ability to select the ideal column for their application. The following section introduces Waters' particle technologies and column formats to help you choose the best column to deliver throughput, resolution, and efficiency for your next chromatographic challenge.

## Particle Technology


Reproducibility and transferability are the cornerstones of Waters' BEH, CSH™, HSS, and solid-core particle technologies. Our extensive portfolio of scalable LC columns exhibit all of the chemical and physical characteristics you would expect from modern LC packing materials.




BEH Technology*	CSH Technology	HSS Technology	Solid-Core Technology
<ul style="list-style-type: none"> <li>High retentivity for basic compounds</li> <li>Exceptional peak shape at elevated pH</li> <li>Good universal column choice for a wide variety of compounds</li> <li>Stable across a wide pH range</li> <li>For separations at high temperatures (80 °C)</li> </ul>	<ul style="list-style-type: none"> <li>Good separations for basic compounds under low pH conditions</li> <li>Excellent MS performance with formic acid as a mobile phase modifier</li> <li>Fast pH switching and column equilibration</li> </ul>	<ul style="list-style-type: none"> <li>High retentivity for polar organic compounds and metabolites</li> <li>Balanced retention of polar and hydrophobic analytes</li> <li>High strength silica for mechanical stability</li> </ul>	<ul style="list-style-type: none"> <li>Maximum efficiency</li> <li>Increased sensitivity</li> <li>Seamless scalability from UPLC to UHPLC to HPLC</li> </ul>

### ETHYLENE BRIDGED HYBRID (BEH) PARTICLE TECHNOLOGY

Ethylene Bridged Hybrid (BEH) columns lead the industry for chromatographic versatility, chemical resistance, and mechanical stability. You can use them at extremes of pH and temperature to enhance retention and specificity for complex mixtures of acidic, basic, and neutral analytes. The BEH-particle family includes general-purpose and application-specific bonded phases that serve application areas.



#### Particle Synthesis




\*US Patents 6,686,035; 7,223,473; 7,250,214.

Refer to "Ethylene-Bridged (BEH Technology™) Hybrids and Their Use in Liquid Chromatography" whitepaper (p/n: [720001159EN](https://www.waters.com/720001159EN)) for further detail.

## CHARGED SURFACE HYBRID (CSH) PARTICLE TECHNOLOGY


Columns packed with charged surface hybrid (CSH) particles manifest the best attributes of BEH particles. CSH stationary phases provide chromatographic selectivity and superior performance in the presence of mobile phases of low ionic strength. The optimized surface charge, pore properties, and bonded phases make charged-surface, hybrid-based columns ideal for rapid method development.

### The Charged-Surface Hybrid Particle



## HIGH STRENGTH SILICA (HSS) PARTICLE TECHNOLOGY

High strength silica (HSS) technology was developed specifically to complement the chromatographic performance of BEH and CSH particles. Compared with the ethylene-bridged BEH and CSH particles, the HSS particle's higher silanophilicity (100% silica) offers chromatographers significant advantages, including increased retention of polar compounds and significantly different selectivity. Additionally, as its name implies, the HSS particle possesses the mechanical strength to operate at pressures as high as 18,000 psi (1240 bar).




## SOLID-CORE PARTICLE TECHNOLOGY

Compared to columns packed with fully porous particles, columns packed with solid-core particles demonstrate higher chromatographic efficiency and lower backpressures. The optimized porous layer that surrounds the solid-silica core gives rise to the key benefits of speed and efficiency. UPLC columns packed with CORTECS™ 1.6 μm particles yield maximum efficiency when used with the ultra-low dispersion ACQUITY UPLC instrument platform. CORTECS Columns packed with 2.7 μm particles offer maximum flexibility, providing increased efficiencies at the backpressure limits of UHPLC and HPLC operation.



## REVERSED PHASE AND HILIC CHEMISTRIES

To maximize selectivity, which is critical for analytical methods development, Waters offers a wide range of chemistries to help you separate the most challenging of compounds. For both reversed-phase and HILIC applicators, Waters has the chemistry you need to get the job done. Rugged and robust ligands that ensure long column lifetimes and reproducible separations over the lifetime of your method. Highlighted are some of Waters chemistries used for reversed-phase and HILIC separations.



## Column Selection

Our quality mission is to ensure that the Waters' Columns you use today are the most reproducible and reliable LC columns available. As a primary manufacturer of silica and hybrid particles, scientists can be assured of consistent column performance, batch-to-batch reproducibility, and product availability over the life of the analytical method.

The following table lists all Waters Column Brands that are registered according to classifications prescribed in the United States Pharmacopeia (USP).

### USP "L" COLUMN LISTING

L1 Octadecyl silane chemically bonded to porous or non-porous silica or ceramic micro-particles, 1.5 to 10 µm in diameter, or monolithic rod			
Brand	Particle Size	Type	Page
AccQTag Ultra RP C <sub>18</sub>	1.7 µm	Spherical	<a href="#">347</a>
ACQUITY UPLC BEH C <sub>18</sub>	1.7 µm	Spherical	<a href="#">118</a> , <a href="#">373</a>
ACQUITY UPLC BEH Shield RP18	1.7 µm	Spherical	<a href="#">118</a>
ACQUITY UPLC CSH C <sub>18</sub>	1.7 µm	Spherical	<a href="#">114</a>
ACQUITY UPLC HSS C <sub>18</sub>	1.7 µm	Spherical	<a href="#">121</a>
ACQUITY UPLC HSS C <sub>18</sub> SB	1.7 µm	Spherical	<a href="#">121</a>
ACQUITY UPLC HSS T3	1.7 µm	Spherical	<a href="#">122</a>
ACQUITY UPLC Oligonucleotide C <sub>18</sub>	1.7 µm	Spherical	<a href="#">377</a>
ACQUITY UPLC Peptide BEH C <sub>18</sub>	1.7 µm	Spherical	<a href="#">377</a>
ACQUITY Premier BEH C <sub>18</sub>	1.7 µm	Spherical	<a href="#">99</a>
ACQUITY Premier HSS T3	187 µm	Spherical	<a href="#">99</a>
ACQUITY Premier Peptide BEH C <sub>18</sub>	1.7 µm	Spherical	<a href="#">100</a>
ACQUITY Premier Oligonucleotide C <sub>18</sub>	1.7 µm	Spherical	<a href="#">100</a>
ACQUITY Premier Shield RP18	1.7 µm	Spherical	<a href="#">99</a>
ACQUITY Premier CSH C <sub>18</sub>	1.7 µm	Spherical	<a href="#">99</a>
XBridge Premier C <sub>18</sub>	2.5 µm	Spherical	<a href="#">101</a>
XBridge Premier Shield RP18	2.5 µm	Spherical	<a href="#">101</a>
XBridge Premier Peptide 130 Å	2.5 µm	Spherical	<a href="#">377</a>
XBridge Premier Peptide 300 Å	2.5 µm	Spherical	<a href="#">377</a>
XBridge Premier Oligonucleotide C <sub>18</sub>	2.5 µm	Spherical	<a href="#">377</a>
XSelect Premier CSH C <sub>18</sub>	2.5 µm	Spherical	<a href="#">101</a>
XSelect Premier HSS T3	2.5 µm	Spherical	<a href="#">102</a>
Atlantis dC <sub>18</sub>	3, 5, 10 µm	Spherical	<a href="#">209</a>
Atlantis T3	3, 5, 10 µm	Spherical	<a href="#">208</a>
CORTECS C <sub>18</sub>	2.7 µm	Spherical	<a href="#">241</a>
CORTECS C <sub>18</sub> +	2.7 µm	Spherical	<a href="#">241</a>
CORTECS Shield RP18	2.7 µm	Spherical	<a href="#">241</a>
CORTECS T3	2.7 µm	Spherical	<a href="#">241</a>
CORTECS UPLC C <sub>18</sub>	1.6 µm	Spherical	<a href="#">110</a>

L1 Octadecyl silane chemically bonded to porous or non-porous silica or ceramic micro-particles, 1.5 to 10 µm in diameter, or a monolithic rod			
Brand	Particle Size	Type	Page
CORTECS UPLC C <sub>18</sub> +	1.6 µm	Spherical	<a href="#">110</a>
CORTECS UPLC Shield RP18	1.6 µm	Spherical	<a href="#">111</a>
CORTECS UPLC T3	1.6 µm	Spherical	<a href="#">111</a>
Delta-Pak C <sub>18</sub>	5 µm	Spherical	<a href="#">301</a>
µBondapak C <sub>18</sub>	10 µm	Irregular	<a href="#">300</a>
µBondapak C <sub>18</sub> Radial-Pak	10 µm	Irregular	<a href="#">307</a>
Nova-Pak C <sub>18</sub>	4, 6 µm	Spherical	<a href="#">299</a>
Prep Nova-Pak HR C <sub>18</sub>	6 µm	Spherical	<a href="#">308</a>
Radial-Pak C <sub>18</sub>	Spherical	Spherical	<a href="#">307</a>
Resolve C <sub>18</sub>	5, 10 µm	Spherical	<a href="#">244</a> , <a href="#">307</a>
Spherisorb ODS1	3, 5, 10 µm	Spherical	<a href="#">295</a>
Spherisorb ODS2	3, 5, 10 µm	Spherical	<a href="#">295</a>
Spherisorb ODS-B	5 µm	Spherical	<a href="#">228</a>
SunFire C <sub>18</sub>	2.5, 3.5, 5, 10 µm	Spherical	<a href="#">172</a> , <a href="#">277</a>
Symmetry C <sub>18</sub>	3.5, 5 µm	Spherical	<a href="#">291</a>
SymmetryPrep C <sub>18</sub>	5, 7 µm	Spherical	<a href="#">217</a> , <a href="#">292</a>
Symmetry 300 C <sub>18</sub>	3.5, 5 µm	Spherical	<a href="#">293</a>
SymmetryShield RP18	3.5, 5 µm	Spherical	<a href="#">292</a> , <a href="#">293</a>
XBridge C <sub>18</sub>	2.5, 3.5, 5, 10 µm	Spherical	<a href="#">255</a>
XBridge Peptide BEH, 130 Å	3.5, 5, 10 µm	Spherical	<a href="#">260</a>
XBridge Peptide BEH, 300 Å	3.5, 5, 10 µm	Spherical	<a href="#">261</a>
XBridge BEH C <sub>18</sub>	2.5, 3.5, 5, 10 µm	Spherical	<a href="#">254</a>
XBridge Oligonucleotide C <sub>18</sub>	2.5 µm	Spherical	<a href="#">262</a>
XBridge Shield RP18	2.5, 3.5, 5, 10 µm	Spherical	<a href="#">256</a>
XSelect CSH C <sub>18</sub>	2.5, 3.5, 5 µm	Spherical	<a href="#">267</a>
XSelect HSS C <sub>18</sub>	2.5, 3.5, 5 µm	Spherical	<a href="#">270</a>
XSelect HSS C <sub>18</sub> SB	2.5, 3.5, 5 µm	Spherical	<a href="#">270</a>
XSelect HSS T3	2.5, 3.5, 5 µm	Spherical	<a href="#">271</a>
XTerra MS C <sub>18</sub>	2.5, 3.5, 5, 10 µm	Spherical	<a href="#">286</a>
XTerra RP18	3.5, 5, 10 µm	Spherical	<a href="#">288</a> , <a href="#">289</a>

( ) - Denotes particle sizes available outside of L class.

Source: United States Pharmacopeia.

L2

Octadecyl silane chemically bonded to silica gel of a controlled surface porosity that has been bonded to a solid spherical core, 30 to 50 µm in diameter

Brand	Particle Size	Type	Page
Bondapak Prep C <sub>18</sub>	15–20 µm	Irregular	<a href="#">236</a>

L3

Porous silica particles, 1.5 to 10 µm in diameter, or a monolithic silica rod

Brand	Particle Size	Type	Page
ACQUITY UPLC BEH HILIC	1.7 µm	Spherical	<a href="#">118,119</a>
Atlantis HILIC Silica	3, 5 µm	Spherical	<a href="#">284</a>
CORTECS HILIC	2.7 µm	Spherical	<a href="#">134</a>
CORTECS UPLC HILIC	1.6 µm	Spherical	<a href="#">112</a>
µPorasil	10 µm	Spherical	<a href="#">237, 244</a>
Nova-Pak Silica	4, 6 µm	Spherical	<a href="#">299</a>
Prep Nova-Pak HR Silica	6 µm	Spherical	<a href="#">308</a>
Resolve Silica	5, 10 µm	Spherical	<a href="#">244, 307</a>
Spherisorb Silica	3, 5, 10 µm	Spherical	<a href="#">230, 295</a>
SunFire Silica	5, 10 µm	Spherical	<a href="#">172, 279</a>
XBridge BEH HILIC	2.5, 3.5, 5, 10 µm	Spherical	<a href="#">265</a>

L4

Silica gel controlled surface porosity bonded to a solid spherical core, 30 to 50 µm in diameter

Brand	Particle Size	Type	Page
Porasil Prep Silica	15–20 µm	Spherical	<a href="#">303</a>

L7

Octylsilane chemically bonded to totally or superficially porous silica particles, 1.5 to 10 µm in diameter, or a monolithic silica rod

Brand	Particle Size	Type	Page
ACQUITY UPLC BEH C <sub>8</sub>	1.7 µm	Spherical	<a href="#">118</a>
CORTECS C <sub>8</sub>	2.7 µm	Spherical	<a href="#">123</a>
CORTECS UPLC C <sub>8</sub>	1.6 µm	Spherical	<a href="#">111</a>
CORTECS Phenyl	2.7 µm	Spherical	<a href="#">112</a>
CORTECS UPLC Phenyl	1.6 µm	Spherical	<a href="#">124</a>
Nova-Pak C <sub>8</sub>	4, 6 µm	Spherical	<a href="#">299</a>
Resolve C <sub>8</sub>	10 µm	Spherical	<a href="#">307</a>
Spherisorb C <sub>8</sub>	3, 5, 10 µm	Spherical	<a href="#">229, 296</a>
SunFire C <sub>8</sub> Silica	3.5, 5, 10 µm	Spherical	<a href="#">173, 280</a>
Symmetry C <sub>8</sub>	3.5, 5, 7 µm	Spherical	<a href="#">292</a>
SymmetryPrep C <sub>8</sub>	7 µm	Spherical	<a href="#">213</a>
SymmetryShield RP8	3.5, 5 µm	Spherical	<a href="#">293</a>
XBridge BEH C <sub>8</sub>	2.5 3.5, 5, 10 µm	Spherical	<a href="#">102, 185</a>
XTerra MS C <sub>8</sub>	2.5 3.5, 5, 10 µm	Spherical	<a href="#">223, 287</a>
XTerra Shield RP8	3.5, 5, 10 µm	Spherical	<a href="#">225, 289</a>

L8

An essentially monomolecular layer of aminopropylsilane chemically bonded to totally porous silica gel support, 1.5 to 10 µm in diameter, or a monolithic silica rod

Brand	Particle Size	Type	Page
High Performance Carbohydrate Analysis	3, 5 µm	–	<a href="#">495</a>
µBondapak NH <sub>2</sub>	10 µm	Irregular	<a href="#">236</a>
Spherisorb NH <sub>2</sub>	3, 5, 10 µm	Spherical	<a href="#">229, 297</a>

L9

Irregular or spherical, totally porous silica gel having a chemically bonded, strongly acidic cation-exchange coating, 3 to 10 µm in diameter

Brand	Particle Size	Type	Page
Spherisorb SCX	5, 10 µm	Spherical	<a href="#">231, 298</a>

L10

Nitrile groups chemically bonded to porous silica particles, 1.5 to 10 µm in diameter, or a monolithic silica rod

Brand	Particle Size	Type	Page
ACQUITY UPLC HSS CN	1.7 µm	Spherical	<a href="#">122</a>
µBondapak CN	10 µm	Irregular	<a href="#">236, 300</a>
NovaPak CN HP	4 µm	Spherical	<a href="#">299</a>
Resolve CN	10 µm	Spherical	<a href="#">307</a>
Spherisorb CN	3, 5, 10 µm	Spherical	<a href="#">230, 297</a>
Spherisorb CN RP	3, 5, 10 µm	Spherical	<a href="#">230, 298</a>
XSelect HSS CN	2.5, 3.5, 5 µm	Spherical	<a href="#">206</a>

L11

Phenyl groups chemically bonded to porous silica particles, 1.5 to 10 µm in diameter, or a monolithic silica rod

Brand	Particle Size	Type	Page
ACQUITY UPLC BEH Phenyl	1.7 µm	Spherical	<a href="#">119</a>
ACQUITY UPLC CSH Phenyl-Hexyl	1.7 µm	Spherical	<a href="#">114</a>
ACQUITY Premier CSH Phenyl-Hexyl	1.7 µm	Spherical	<a href="#">99</a>
XSelect Premier CSH Phenyl-Hexyl	1.7 µm	Spherical	<a href="#">102</a>
CORTECS Phenyl	2.7 µm	Spherical	<a href="#">136</a>
CORTECS UPLC Phenyl	1.6 µm	Spherical	<a href="#">112</a>
µBondapak Phenyl	10 µm	Irregular	<a href="#">236</a>
NovaPak Phenyl	4 µm	Spherical	<a href="#">233</a>
Spherisorb Phenyl	3, 5, 10 µm	Spherical	<a href="#">230</a>
XBridge BEH Phenyl	2.5, 3.5, 5 µm	Spherical	<a href="#">147</a>
XSelect CSH Phenyl-Hexyl	2.5, 3.5, 5 µm	Spherical	<a href="#">158</a>
XTerra Phenyl	3.5, 5 µm	Spherical	<a href="#">177</a>

L12

A strong anion-exchange packing made by chemically bonding a quaternary amine to a solid silica spherical core, 30 to 50 µm in diameter

Brand	Particle Size	Type	Page
AccellPlus QMA	40 µm	Irregular	<a href="#">450</a>

( ) - Denotes particle sizes available outside of L class.

Source: United States Pharmacopeia.

<b>L13</b> Trimethylsilane chemically bonded to porous silica particles, 3 to 10 µm in diameter			
Brand	Particle Size	Type	Page
Spherisorb C <sub>1</sub>	3, 5, 10 µm	Spherical	<a href="#">229</a> , <a href="#">296</a>

<b>L14</b> Silica gel having a chemically bonded strongly basic quaternary ammonium anion-exchange coating, 5 to 10 µm in diameter			
Brand	Particle Size	Type	Page
Spherisorb SAX	5, 10 µm	Spherical	<a href="#">231</a> , <a href="#">298</a>

<b>L15</b> Hexylsilane chemically bonded to totally porous silica particles, 3 to 10 µm in diameter			
Brand	Particle Size	Type	Page
Spherisorb C <sub>6</sub>	3, 5, 10 µm	Spherical	<a href="#">229</a> , <a href="#">296</a>

<b>L17</b> Strong cation-exchange resin consisting of sulfonated cross-linked styrene-divinylbenzene copolymer in the hydrogen form, 6 to 12 µm in diameter			
Brand	Particle Size	Type	Page
Fast Fruit Juice	N/A	N/A	<a href="#">497</a>
IC-Pak Cation	10 µm	Irregular	<a href="#">500</a>
IC-Pak Ion Exclusion	7 µm	Spherical	<a href="#">497</a>
Shodex RSPak DC-613	6 µm	Spherical	<a href="#">237</a>

<b>L19</b> Strong cation-exchange resin consisting of sulfonated cross-linked styrene-divinylbenzene copolymer in the calcium form, 5 to 15 µm in diameter			
Brand	Particle Size	Type	Page
Shodex Sugar SC-1011	7 µm	Spherical	<a href="#">496</a>
Sugar-Pak 1	9 µm	Spherical	<a href="#">496</a>

<b>L20</b> Dihydroxypropane groups chemically bonded to porous silica or hybrid particles, 1.5 to 10 µm in diameter, or a monolithic silica rod			
Brand	Particle Size	Type	Page
ACQUITY BEH200SEC	1.7 µm	Spherical	<a href="#">119</a>
BioSuite 125, 250, 450 series	4, 5, 8, 10, (13), (17) µm	Spherical	<a href="#">433</a>
Insulin HMWP	-	N/A	<a href="#">429</a>
Protein-Pak 60	10 µm	Spherical	<a href="#">433</a>
Protein-Pak 125	10 µm	Spherical	<a href="#">433</a>
Protein-Pak 200SW and 300SW	10 µm	Spherical	<a href="#">433</a>
XBridge Protein BEH SEC, 125 Å	3.5 µm	Spherical	<a href="#">466</a>
XBridge Protein BEH SEC, 200 Å	3.5 µm	Spherical	<a href="#">466</a>
XBridge Protein BEH SEC, 450 Å	3.5 µm	Spherical	<a href="#">466</a>

<b>L21</b> A rigid, spherical styrene-divinylbenzene copolymer, 3 to 30 µm in diameter			
Brand	Particle Size	Type	Page
Styragel HR 0.5, 1, 2, 3 and 4	-	Spherical	<a href="#">406</a>
Styragel HR 4E	-	Spherical	<a href="#">406</a>
Styragel HR 5E	-	Spherical	<a href="#">406</a>

<b>L22</b> A cation-exchange resin made of porous polystyrene gel with sulfonic acid groups, 5 to 15 µm in diameter			
Brand	Particle Size	Type	Page
IC-Pak Ion Exclusion	7 µm	Spherical	<a href="#">500</a>
Shodex RSPak DC-613	6 µm	Spherical	<a href="#">237</a>
Shodex Sugar SP0810	8 µm	Spherical	<a href="#">496</a>

<b>L23</b> An anion-exchange resin made of porous polymethacrylate or polyacrylate gel with quaternary ammonium groups, 7 to 12 µm in size			
Brand	Particle Size	Type	Page
BioSuite DEAE	(2.5), 10, 13 µm	Spherical	<a href="#">446</a>
BioSuite Q AXC	10, 13 µm	Spherical	<a href="#">446</a>
BioSuite Q-PEEK	10 µm	Spherical	<a href="#">446</a>
IC-Pak Anion	10 µm	Spherical	<a href="#">500</a>
IC-Pak A HC	10 µm	Spherical	<a href="#">500</a>
Protein-Pak Q 8HR	8 µm	Spherical	<a href="#">447</a>

<b>L25</b> Packing having the capacity to separate compounds with a molecular weight range from 100-5000 (as determined by polyethylene oxide), applied to neutral, anionic, and cationic water-soluble polymers. A polymethacrylate resin base, cross-linked with polyhydroxylated ether (surface contained some residual carboxyl functional groups) was found suitable			
Brand	Particle Size	Type	Page
Ultrahydrogel DP, +120	10 µm	Spherical	<a href="#">473</a>

<b>L26</b> Butyl silane chemically bonded to totally porous or superficially porous silica particles, 1.5 to 10 µm in diameter			
Brand	Particle Size	Type	Page
ACQUITY UPLC BEH300 C <sub>4</sub>	1.7 µm	Spherical	<a href="#">416</a>
Delta-Pak C <sub>4</sub>	5 µm	Spherical	<a href="#">234</a>
Symmetry300 C <sub>4</sub>	3.5, 5 µm	Spherical	<a href="#">241</a> , <a href="#">242</a>
XBridge BEH300 C <sub>4</sub>	3.5, 5, 10 µm	Spherical	<a href="#">241</a>

<b>L27</b> Porous silica particles, 30 to 50 µm in diameter			
Brand	Particle Size	Type	Page
Porasil	37-55 µm	Spherical	<a href="#">303</a> , <a href="#">307</a>

<b>L33</b> Packing having the capacity to separate dextrans by molecular size over a range of 4000 to 500,000 Da. It is spherical, silica-based, and processed to provide pH stability			
Brand	Particle Size	Type	Page
ACQUITY UPLC Protein BEH SEC, 125 Å	1.7 µm	Spherical	<a href="#">427</a>

(-) - Denotes particle sizes available outside of L class.

Source: United States Pharmacopeia.

**L34** Strong cation-exchange resin consisting of sulfonated cross-linked styrene-divinylbenzene copolymer in the lead form, 7 to 9 µm in diameter

Brand	Particle Size	Type	Page
Shodex Sugar SP-0810	N/A	Spherical	<a href="#">496</a>

**L37** Packing having the capacity to separate proteins by molecular size over a range of 2000 to 40,000 Da. It is a polymethacrylate gel

Brand	Particle Size	Type	Page
Ultrahydrogel 250	N/A	Spherical	<a href="#">473</a>

**L38** A methacrylate-based size-exclusion packing for water-soluble samples

Brand	Particle Size	Type	Page
Ultrahydrogel series	N/A	Spherical	<a href="#">473</a>

**L39** A hydrophilic polyhydroxymethacrylate gel of totally porous spherical resin

Brand	Particle Size	Type	Page
Ultrahydrogel series	N/A	Spherical	<a href="#">473</a>

**L43** Pentafluorophenyl groups chemically bonded to silica particles by a propyl spacer, 1.5 to 10 µm in diameter

Brand	Particle Size	Type	Page
ACQUITY UPLC CSH Fluoro-Phenyl	1.7 µm	Spherical	<a href="#">114</a>
ACQUITY UPLC HSS PFP	1.8 µm	Spherical	<a href="#">122</a>
XSelect CSH Fluoro-Phenyl	2.5, 3.5, 5 µm	Spherical	<a href="#">157, 160, 161</a>
XSelect HSS PFP	2.5, 3.5, 5 µm	Spherical	<a href="#">165, 167</a>

**L52** A strong cation exchange resin made of porous silica with sulfopropyl or sulfoethyl groups, 1 to 10 µm in diameter

Brand	Particle Size	Type	Page
IC-Pak Cation	10 µm	Irregular	<a href="#">500</a>

**L55** A strong cation exchange resin made of porous silica coated with polybutadiene-maleic acid copolymer, about 5 µm in diameter

Brand	Particle Size	Type	Page
IC-Pak C M/D	N/A	N/A	<a href="#">500</a>

**L59** Packing for the size-exclusion separations of proteins (separation by molecular weight) over the range of 5 to 7000 kDa. The packing is spherical 1.5 to 10 µm, silica or hybrid packing with a hydrophilic coating

Brand	Particle Size	Type	Page
ACQUITY BEH200 SEC	1.7 µm	Spherical	<a href="#">119</a>
BioSuite 125, 250, 450 series	4–17 µm	Spherical	<a href="#">433</a>
Protein-Pak 60	10 µm	Spherical	<a href="#">433, 471</a>
Protein-Pak 300SW	10 µm	Spherical	<a href="#">433, 471</a>

**L68** Spherical, porous silica, 10 µm or less in diameter, the surface of which has been covalently modified with alkyl amide groups and not endcapped

Brand	Particle Size	Type	Page
ACQUITY UPLC Glycan BEH Amide	1.7 µm	Spherical	<a href="#">119</a>
ACQUITY UPLC BEH Amide	1.7 µm	Spherical	<a href="#">118</a>
ACQUITY Premier Glycan BEH Amide	1.7 µm	Spherical	<a href="#">369</a>
ACQUITY Premier BEH Amide	1.7 µm	Spherical	<a href="#">383</a>
XBridge Premier Glycan BEH Amide	2.5 µm	Spherical	<a href="#">338</a>
XBridge BEH Amide	2.5, 3.5, 5 µm	Spherical	<a href="#">338</a>
XBridge Premier BEH Amide	2.5 µm	Spherical	<a href="#">101</a>
XBridge BEH Amide Glycan	2.5, 3.5 µm	Spherical	<a href="#">370</a>

**L78** A silane ligand that consists of both reversed-phase (an alkyl chain longer than C<sub>8</sub>) and anion-exchange (primary, secondary, tertiary, or quaternary amino groups) functional groups chemically bonded to porous or non-porous silica or ceramic micro-particles, 1.0 to 50 µm in diameter, or a monolithic rod.

Brand	Particle Size	Type	Page
Atlantis Premier BEH C <sub>18</sub> AX	1.7, 2.5, 5 µm	Spherical	<a href="#">104</a>

**L122** Sulfobetaine graft-polymerized to totally or superficially porous hydrophilic polymer particles, 1.0 to 10 µm in diameter, or a monolithic rod. Packing having densely bonded zwitterionic groups with 1:1 charge balance

Brand	Particle Size	Type	Page
Atlantis Premier BEH Z-HILIC	1.7, 2.5, 5 µm	Spherical	<a href="#">104</a>

( ) - Denotes particle sizes available outside of L class.

Source: United States Pharmacopeia.



## Column Configurations for Any LC System

### COLUMN NOMENCLATURE

Our fully scalable particle technologies ensure that our LC columns perform with a broad range of chromatographic instrumentation. Depending on the goals of a separation, the instrument platform used, or the sample type, you can choose the most suitable column that is matched to your system's configuration without adversely affecting the chromatographic result.

The following table serves as a guide for selecting an appropriate LC column according to instrument classification.

Nano/Micro	UPLC	UHPLC	HPLC	Preparative
ACQUITY UPLC M-CLASS BEH (1.7 µm)	ACQUITY UPLC BEH (1.7 µm)	XBridge BEH <i>XP</i> (2.5 µm)	XBridge BEH (3.5, 5 µm)	XBridge BEH OBD™ (5, 10 µm)
ACQUITY UPLC M-CLASS CSH (1.7 µm)	ACQUITY UPLC CSH (1.7 µm)	XSelect CSH <i>XP</i> (2.5 µm)	XSelect CSH <i>XP</i> (3.5, 5 µm)	XSelect CSH OBD (5, 10 µm)
ACQUITY UPLC M-CLASS HSS (1.8 µm)	ACQUITY UPLC HSS (1.8 µm)	XSelect HSS <i>XP</i> (2.5 µm)	XSelect HSS <i>XP</i> (3.5, 5 µm)	XSelect HSS OBD (5 µm)
—	CORTECS UPLC (1.6 µm)	CORTECS (2.7 µm)	—	—

### COLUMN CONFIGURATION

System dispersion is inherent in every chromatographic system. It is the instrument's contribution to chromatographic band broadening and is dependent on the system's tubing volume, valve fittings, column fittings, and flow cell volume. System dispersion, in combination with column dispersion, makes up the total dispersion of a given separation. Therefore, it is important to understand the system's impact on chromatographic band broadening when choosing your column configuration. Systems that have high dispersion values will obtain the best column performance using columns that have larger volumes; and, systems that have low dispersion values are able to obtain excellent column performance using columns that have smaller volumes.

The following table summarizes the characteristics of Waters LC Systems and matches the column configuration that maintains chromatographic efficiency.



System	Nano/Micro	UPLC	UHPLC	HPLC	Preparative
Dispersion	1 µL	<20 µL	22–29 µL	>30 µL	—
Routine Pressure	<15,000 psi	<18,000 psi	<10,000 psi	<10,000 psi	<4000 psi
Particle Size	<2 µm	<2 µm	2–3 µm	3–5 µm	>5 µm
Column I.D.	75–300 µm	2.1 mm (1.0 mm)	3.0 mm (2.1 mm)	4.6 mm (3.0 mm)	>7.8 mm
Column Length	50–250 mm	<150 mm	50–150 mm	75–300 mm	50–300 mm

When you transfer LC methods, instrument dispersion is one of the most practical LC-instrument parameters to determine. Knowing the bandsread value helps you develop your own compatible methods, allowing you to seamlessly scale column dimensions or transfer methods between different instrumentation platforms and laboratory functions. The following table recommends column configurations based on nominal instrument bandsread values.

System	Bandsread*	Recommended Column Particle Sizes and I.D.s	System	Bandsread*	Recommended Column Particle Sizes and I.D.s
Shimadzu Prominence UFLC	41 µL	CORTECS 2.7 µm	ACQUITY UPLC	12 µL	ACQUITY UPLC BEH 1.7 µm
Alliance 2695 HPLC	29 µL	XBridge 3.5, 5 µm	ACQUITY UPLC H-Class with Column Manager	12 µL	ACQUITY UPLC CSH 1.7 µm
Agilent 1260 UHPLC (600 bar)	28 µL	XSelect 3.5, 5 µm			ACQUITY UPLC HSS 1.8 µm
		<b>3.0–4.6 mm I.D.</b>	ACQUITY UPLC H-Class	9 µL	CORTECS UPLC 1.6 µm
ACQUITY Arc	23 µL	XBridge 2.5, 3.5, 5 µm			<b>2.1 mm I.D.</b>
		XSelect 2.5, 3.5, 5 µm	ACQUITY UPLC I-Class (FTN)	7.5 µL	ACQUITY UPLC BEH 1.7 µm
		CORTECS 2.7 µm			ACQUITY UPLC CSH 1.7 µm
		<b>3.0 mm I.D.</b>	ACQUITY UPLC I-Class (FL)	5.5 µL	ACQUITY UPLC HSS 1.8 µm
Thermo Accela UHPLC	21 µL	XBridge 2.5, 3.5, 5 µm			CORTECS UPLC 1.6 µm
		XSelect 2.5, 3.5, 5 µm			<b>1.0–2.1 mm I.D.</b>
Agilent 1290 UHPLC (1200 bar)	17 µL	CORTECS 2.7 µm			
		<b>3.0 mm I.D.</b>			

\*These data are based on nominal values for unmodified systems. As such, they are intended for reference only. Any adjustment to a system's plumbing, connectivity, and configuration will change the instrument's bandsread, affecting the quality of chromatography.

### L/d<sub>p</sub> COMPARISON CHART FOR LC COLUMNS

To convert an HPLC method to a UPLC or UHPLC method with no loss in resolution, select columns that have equivalent length-to-particle-size (L/d<sub>p</sub>) ratio.

Waters uses this ratio to compare the resolving power of columns. If you keep the L/d<sub>p</sub> ratio the same for two columns, you will obtain the same resolution. Therefore, for two columns with the same L/d<sub>p</sub> ratio, the more efficient, shorter column (packed with smaller particles) will provide the same resolution in less time.

$$\text{EXAMPLE: } \frac{150 \text{ mm}}{5 \text{ }\mu\text{m}} = \frac{150,000 \text{ }\mu\text{m}}{5 \text{ }\mu\text{m}} = 30,000$$

L/d <sub>p</sub>		Column length (mm)							
		20	30	50	75	100	150	250	
Particle size (µm)	Fully porous	1.7	-	17,600	29,400	44,100	58,000	88,200	-
		2.5	8,000	12,000	20,000	30,000	40,000	60,000	-
		3.5	5,700	8,600	14,300	21,400	28,600	42,900	71,400
		5.0	4,000	6,000	10,000	15,000	20,000	30,000	50,000