




CHIRAL
TECHNOLOGIES INC


SUBSIDIARY OF DAICEL CHEMICAL INDUSTRIES, LTD.

CHIRALPAK® AY™

***New proprietary
polysaccharide chiral
stationary phases and columns —
another perfect solution.***







- *Authentic Daicel stationary phases and columns*
- *Unique selectivity*
- *Proven Daicel quality*
- *Excellent resolution power*
- *Wide variety of column dimensions and particle sizes*



AUTHORIZED DISTRIBUTOR

MZ-Analysentechnik GmbH
Barcelona-Allee 17 • D-55129 Mainz
Tel +49 6131 880 96-0
Fax +49 6131 880 96-20
e-mail: info@mz-at.de
www.mz-at.de

CHIRALPAK® AY™



For more information about all of our proprietary columns and CSPs, call 1-800-6CHIRAL, or visit us online at www.chiraltech.com.

PART NO.	PRODUCT NAME	PARTICLE SIZE (µ)	ID (mm)	LENGTH (mm)	PRODUCT TYPE
47020	CHIRALPAK AY	20	1 kg	–	CSP
47225	CHIRALPAK AY	20	4.6	250	Method Develop
47233	CHIRALPAK AY	20	10	100	SMB (set of 8)
47242	CHIRALPAK AY	20	21	50	Semi-Prep Guard
47245	CHIRALPAK AY	20	21	250	Semi-Prep
47256	CHIRALPAK AY	20	50	500	Preparative
47266	CHIRALPAK AY	20	100	500	Preparative

PART NO.	PRODUCT NAME	PARTICLE SIZE (µ)	ID (mm)	LENGTH (mm)	PRODUCT TYPE
473C4	CHIRALPAK AY-H	5	0.3	150	Microflow
47311	CHIRALPAK AY-H	5	4.0	10	Guard Cartridges
47324	CHIRALPAK AY-H	5	4.6	150	Analytical
47325	CHIRALPAK AY-H	5	4.6	250	Analytical
47394	CHIRALPAK AY-H	5	2.1	150	Analytical