




CHIRAL
TECHNOLOGIES INC

SUBSIDIARY OF  DAICEL CHEMICAL INDUSTRIES, LTD.


CHIRALPAK® AZ™

***New proprietary
polysaccharide chiral
stationary phases and columns —
the perfect solution.***

R=







- *Authentic Daicel stationary phases and columns*
- *Unique selectivity*
- *Proven Daicel quality*
- *Excellent resolution power*
- *Wide variety of column dimensions and particle sizes*



AUTHORIZED DISTRIBUTOR

MZ-Analysentechnik GmbH
Barcelona-Allee 17 • D-55129 Mainz
Tel +49 6131 880 96-0
Fax +49 6131 880 96-20
e-mail: info@mz-at.de
www.mz-at.de



For more information about all of our proprietary columns and CSPs, call 1-800-6CHIRAL, or visit us online at www.chiraltech.com.

PART NO.	PRODUCT NAME	PARTICLE SIZE (µ)	ID (mm)	LENGTH (mm)	PRODUCT TYPE
61161	CHIRALPAK AZ	20	1 kg	—	CSP
61225	CHIRALPAK AZ	20	4.6	250	Method Develop
61233	CHIRALPAK AZ	20	10	100	SMB (set of 8)
61235	CHIRALPAK AZ	20	10	250	Semi-Prep
61242	CHIRALPAK AZ	20	21	50	Semi-Prep Guard
61245	CHIRALPAK AZ	20	21	250	Semi-Prep
61255	CHIRALPAK AZ	20	50	250	Preparative
61256	CHIRALPAK AZ	20	50	500	Preparative
61266	CHIRALPAK AZ	20	100	500	Preparative
61275	CHIRALPAK AZ	20	30	250	Preparative

PART NO.	PRODUCT NAME	PARTICLE SIZE (µ)	ID (mm)	LENGTH (mm)	PRODUCT TYPE
61935	CHIRALPAK AZ	20	10	250	SFC Semi-Prep
61945	CHIRALPAK AZ	20	20	250	SFC Semi-Prep
61955	CHIRALPAK AZ	20	50	250	SFC Semi-Prep
61975	CHIRALPAK AZ	20	30	250	SFC Semi-Prep
613C4	CHIRALPAK AZ-H	5	0.3	150	Microflow
61311	CHIRALPAK AZ-H	5	4.0	10	Guard Cartridge
61324	CHIRALPAK AZ-H	5	4.6	150	Analytical
61325	CHIRALPAK AZ-H	5	4.6	250	Analytical
61423	CHIRALPAK AZ-H	5	4.6	100	SFC Screening
61724	CHIRALPAK AZ-RH	5	4.6	150	Analytical