




CHIRAL
TECHNOLOGIES INC

SUBSIDIARY OF  DAICEL CHEMICAL INDUSTRIES, LTD.

CHIRALPAK® IA™

***Immobilized Polysaccharide
chiral stationary phases
and columns — the solution
you've been looking for.***





- *Authentic Daicel stationary phases and columns*
- *Outstanding solvent tolerance*
- *Solvent flexibility for the resolution of compounds with limited solubility*
- *Extended durability*
- *Excellent separation of enantiomers*



AUTHORIZED DISTRIBUTOR

MZ-Analysentechnik GmbH
Barcelona-Allee 17 • D-55129 Mainz
Tel +49 6131 880 96-0
Fax +49 6131 880 96-20
e-mail: info@mz-at.de
www.mz-at.de

CHIRALPAK® IA™



All UV at 254 nm

For more information about all of our proprietary columns and CSPs, call 1-800-6CHIRAL, or visit us online at www.chiraltech.com.

PART NO.	PRODUCT NAME	PARTICLE SIZE (µ)	ID (mm)	LENGTH (mm)	PRODUCT TYPE
80020	CHIRALPAK IA	20	1 kg	—	CSP
80225	CHIRALPAK IA	20	4.6	250	Method Develop
80233	CHIRALPAK IA	20	10	100	SMB (set of 8)
80242	CHIRALPAK IA	20	21	50	Semi-Prep Guard
80245	CHIRALPAK IA	20	21	250	Semi-Prep
80256	CHIRALPAK IA	20	50	500	Preparative
80266	CHIRALPAK IA	20	100	500	Preparative
803C4	CHIRALPAK IA	5	0.3	150	Microflow
80311	CHIRALPAK IA	5	4.0	10	Guard Cartridges
80324	CHIRALPAK IA	5	4.6	150	Analytical
80325	CHIRALPAK IA	5	4.6	250	Analytical
80335	CHIRALPAK IA	5	10	250	Semi-Prep

PART NO.	PRODUCT NAME	PARTICLE SIZE (µ)	ID (mm)	LENGTH (mm)	PRODUCT TYPE
80345	CHIRALPAK IA	5	20	250	Semi-Prep
80375	CHIRALPAK IA	5	30	250	Semi-Prep
80394	CHIRALPAK IA	5	2.1	150	Analytical
80422	CHIRALPAK IA	5	4.6	50	SFC Guard
80423	CHIRALPAK IA	5	4.6	100	SFC Screening
80435	CHIRALPAK IA	5	10	250	SFC Semi-Prep
80445	CHIRALPAK IA	5	21	250	SFC Semi-Prep
80455	CHIRALPAK IA	5	50	250	SFC Preparative
80475	CHIRALPAK IA	5	30	250	SFC Semi-Prep
80522	CHIRALPAK IA-3	3	4.6	50	Analytical
80523	CHIRALPAK IA-3	3	4.6	100	Analytical
80524	CHIRALPAK IA-3	3	4.6	150	Analytical